[Gebrauchsanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) |
[Instrucciones de manejo](#) | [Istruzione](#) | [Instruções de utilização](#) | [操作手册](#) |
[Руководство по эксплуатации](#) | [사용 지침](#) | [Bruksanvisning](#) | [使用説明書](#) |
[Használati utasítás](#) | [Návod k použití](#) | [Gebruiksaanwijzing](#) |
[Instrukcja użytkowania](#) | [Kullanım Talimatları](#)

Dispensette® S

Dispensette® S Organic

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Table of Contents

1 Introduction	52	15 Repairs	95
1.1 Scope of delivery	52	15.1 Sending for repair	95
1.2 Use operating manual	52	16 Calibration service	96
2 Safety regulations.....	54	17 Warranty.....	97
2.1 General safety regulations	54	18 Disposal	98
2.2 Target group.....	55		
2.3 Improper use	55		
2.4 Foreseeable misuse	55		
2.5 Function.....	55		
2.6 Limitations of Use	56		
2.7 Operating Limitations.....	57		
2.8 Operating exclusions	57		
2.9 Storage conditions.....	57		
2.10 Dispenser Selection Chart	58		
3 Functional and operational components	60		
4 Assembly.....	62		
4.1 First Steps.....	62		
4.2 Priming	63		
5 Operation.....	65		
5.1 Dispensing.....	65		
5.2 Accessories	66		
6 Error limits	69		
7 Checking the Volume (Calibration)	70		
8 Adjustment.....	71		
8.1 Digital model	71		
8.2 Analog-adjustable model	72		
9 Cleaning	74		
9.1 Cleaning type Analog, type Fixed	74		
9.2 Cleaning Digital model	76		
9.3 Replacing the dosing cannula/ valves.....	79		
10 Autoclaving	83		
11 Troubleshooting	84		
12 Product markings	86		
13 Ordering Information	87		
14 Accessories/spare parts	89		

1 Introduction

1.1 Scope of delivery

Bottle-top dispenser Dispensette® S or Dispensette® S Organic, for GL 45 thread bottles, discharge tube or discharge tube with recirculation valve, telescopic filling tube (optional for instruments with recirculation valve), installation tool, various bottle adapters, a performance certificate and these operating instructions.

Nominal volume ml	Adapter for bottle thread, PP	Filling tube length mm
1, 2, 5, 10	GL 24-25, GL 28/S 28, GL 32-33, GL 38, S 40	125–240
25, 50, 100	GL 32-33, GL 38, S 40	170–330

1.2 Use operating manual

- Please carefully read the operating manual before using the device for the first time.
- Keep the operating manual in an easily accessible place. It is part of the instrument.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find updated versions of the Instructions for Use at www.brand.de/om.

1.2.1 Signal words and their meaning

Signal word	Meaning
 WARNING or  WARNING! ...	WARNING indicates a dangerous situation that, if not avoided, could result in death or serious injury.
 CAUTION or  CAUTION! ...	CAUTION indicates a hazardous situation that, if not avoided, could result in moderate or minor injury.
NOTICE or NOTICE ...	NOTE is used to address actions that are not related to physical injury. Example: Possible property damage.

1.2.2 Symbols used in the instructions

The sign warns of a possible hazard.

The sign suggests that eye protection be worn.

The sign warns of hazardous media and the dangers that can arise from liquid splashes.

The sign suggests that you to wear protective clothing.

The symbol indicates that the instruction manual and safety information must be observed.

The sign suggests that protective gloves be worn.

The sign suggests that a face shield be worn during work.

1.2.3 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Safety regulations

2.1 General safety regulations

Please read carefully!

The Dispensette® S can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

- Every user must read and observe the user manual enclosed with this laboratory device before and during use. The device may only be used by trained and qualified personnel.
- Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves.
- When working with infectious or hazardous samples/media (e.g. hazardous substances), the general safety rules in the laboratory must be observed and regulations for handling the samples/media must be followed. The information provided by the media manufacturer (e.g. safety data sheets) must be observed.
- Use the laboratory device only for pipetting or dispensing media within the defined limitations and restrictions of use. Comply with the operating exclusions.
- If working with inflammable media, make sure to take the necessary precautions to avoid the buildup of static electricity (e.g., do not dispense into plastic vessels and do not wipe instruments with a dry cloth). Do not use the device in potentially explosive atmospheres. If in doubt, contact the manufacturer or supplier.
- Always check that the instrument is in proper working condition before use. If device malfunctions are indicated (e.g. sluggish pistons, leaks or power supply faults), stop working with the device immediately and refer to the "Troubleshooting" section in the user manual. Contact the manufacturer, if necessary.
- Always perform work in a manner that does not endanger users or other people. Avoid splashes. Only use suitable vessels. Never use unnecessary force when operating, cleaning or servicing the laboratory device.
- If the laboratory device is supplied with power through a power adapter, batteries or accumulators, the components and connections on the device must be regularly checked for proper working condition. Do not operate the laboratory device and its accessories in an unprotected, humid or wet environment.
- Do not make any technical modifications. Use only original replacement parts of the manufacturer; even power adapters or accumulators of identical size and specification from other manufacturers are not permitted. Do not disassemble the laboratory device and its accessories (e.g., power adapters, cables, racks, accumulators or batteries) any further than described in the user manual!
- Only autoclave the laboratory device when this is permissible according to the user manual.
- Never press down the piston when the dispensing cannula is sealed with the screw cap.
- Never remove the dispensing cannula while the dispensing cylinder is filled.
- Reagents can accumulate in the screw cap of the dispensing cannula. Therefore, clean the screw cap regularly.
- For small bottles, and when using the flexible discharge tube, use a bottle stand to prevent tipping over.
- Never carry the mounted instrument by the cylinder sleeve or the valve block. Breakage and detachment of the cylinder can lead to injuries caused by chemicals (First Steps, p. 62, Fig. 3).
- Never use force. Always pull the piston up and press down gently when dispensing.

2.2 Target group

The operating manual is intended for users who use the laboratory instrument in the course of their professional activities. Users are familiar with the typical safety regulations and working methods in laboratories and have been trained accordingly. They can recognize potential hazards and protect themselves from them. The operating manual assumes this expertise and does not replace basic laboratory training or specific safety training.

2.3 Improper use

Various risks may arise if the laboratory instrument is used improperly. These risks include: inaccurate liquid delivery, damage to the laboratory instrument, and the risk of contamination, infection, and injury from contact with the dispensed media.

Any use of the laboratory instrument other than for dispensing liquids within the defined operating limits is considered improper use.

2.4 Foreseeable misuse

A typical misuse is the incomplete dispensing of liquids, i.e. residual liquid remains in the instrument if the piston is not pressed down completely. If this is the case, liquid will escape from the upper edge of the dispensing cylinder over time due to the design.

Another typical misuse is inadequate cleaning of crystallizing or hygroscopic liquids. This leads to blocked pistons or valves. In the case of hygroscopic liquids, this leads to an accumulation of liquid at the recirculation valve and in the dispensing cylinder.

2.5 Function

With the Dispensette® S and Dispensette® S Organic bottle-top dispensers, liquids can be dispensed directly from the supply bottle. Available in Digital, Analog and Fixed models. The instruments are marked DE-M and optionally equipped with recirculation valve.


- Dispensette® S offers a wide range of applications for dispensing aggressive reagents (e.g., concentrated acids such as H_3PO_4 and H_2SO_4 , bases such as NaOH and KOH, and salt solutions as well as various polar solvents).
- Dispensette® S Organic is used for dispensing organic solvents (e.g., chlorinated and fluorinated hydrocarbons such as trichlorotrifluoroethane and dichloromethane or acids such as concentrated HCl and HNO_3 as well as for trifluoroacetic acid (TFA), tetrahydrofuran (THF) and peroxides).

NOTICE For guidelines on selecting the right dispenser, please observe the corresponding exclusions of use and the following “Dispenser Selection Chart”.

NOTICE **Dispensing hydrofluoric acid (HF):** For dispensing hydrofluoric acid (HF), we recommend the bottle-top dispenser Dispensette® S Trace Analysis with platinum-iridium valve spring (see separate operating instructions).

2.5.1 Color codes

Dispensette® S

2.5.2 Handling

If handled properly, the dispensed liquid comes into contact with the following chemical-resistant materials only:

Dispensette® S	Borosilicate glass, Al ₂ O ₃ -ceramic, ETFE, FEP, PFA, PTFE, platinum-iridium, PP (screw cap, bottle adapter).
Dispensette® S Organic	Borosilicate glass, Al ₂ O ₃ -ceramic, ETFE, FEP, PFA, PTFE, tantalum, PP (screw cap, bottle adapter).

NOTICE Bottle adapter: Alternatively, ETFE/PTFE bottle adapters can be used (Accessories/spare parts, p. 89). The user must review carefully to ensure the suitability of ETFE/PTFE bottle adapters.

NOTICE Dispensing hydrofluoric acid (HF): For dispensing hydrofluoric acid (HF), we recommend the bottle-top dispenser Dispensette® S Trace Analysis with platinum-iridium valve spring (see separate operating instructions).

2.6 Limitations of Use

This instrument is designed for dispensing liquids, observing the following physical limits:

- operating temperature from +15°C to +40°C (from 59°F to 104°F) of instrument and reagent
- vapor pressure up to 600 mbar. Absorb slowly above 300 mbar in order to prevent the liquid from boiling
- Kinematic viscosity up to 500 mm²/s (dynamic viscosity [mPas] = kinematic viscosity [mm²/s] x density [g/cm³])
- Density up to 2.2 g/cm³

2.7 Operating Limitations

- Liquids that form deposits (e.g., crystallizing solutions or concentrated alkaline solutions) may make the piston difficult to move or may cause jamming. If the piston movement becomes sluggish or stiff, the instrument should be cleaned immediately. See also Cleaning, p. 74.
- When dispensing inflammable media, make sure to avoid to buildup of static charge (e.g., do not dispense into plastic vessels; do not wipe instruments with a dry cloth).
- The instrument is designed for general laboratory applications and complies with the requirements of the relevant standards, e.g. DIN EN ISO 8655. The use of the instrument for special applications (e.g. in trace analysis, in the food sector, etc.) must be carefully checked by the user. Approvals for specific applications (e.g., for the production and administration of food, pharmaceuticals, or cosmetics) have not been given.

2.8 Operating exclusions

Never use Dispensette® S for:

- Liquids that attack Al_2O_3 -ceramic, ETFE, FEP, PFA, and PTFE (e.g., dissolved sodium azide*)
- Liquids that attack borosilicate glass (e.g., hydrofluoric acid)
- Liquids that are decomposed catalytically by platinum-iridium (e.g., H_2O_2)
- Hydrochloric acid > 20% and nitric acid > 30%
- Tetrahydrofuran
- Trifluoroacetic acid
- Explosive liquids
- Carbon disulfide
- Suspensions (e.g., activated charcoal) because solid particles may clog or damage the instrument
- Liquids that tend to polymerize spontaneously (e.g., adhesives and their components)
- Liquids that attack PP (screw cap and bottle adapter)**.

Never use Dispensette® S Organic for:

- Liquids that attack Al_2O_3 -ceramic, tantalum, ETFE, FEP, PFA, and PTFE (e.g., dissolved sodium azide*)
- Liquids that attack borosilicate glass (e.g., hydrofluoric acid)
- Bases and salt solutions
- Explosive liquids
- Carbon disulfide
- Suspensions (e.g., activated charcoal) because solid particles may clog or damage the instrument
- Liquids that tend to polymerize spontaneously (e.g., adhesives and their components)
- Liquids that attack PP (screw cap and bottle adapter)**.

* Dissolved sodium azide is permitted up to a concentration of 0.1%.

** Alternatively, ETFE/PTFE bottle adapters and ETFE/PTFE screw caps can be used (accessories from Accessories/spare parts, p. 89). The user must review carefully to ensure the suitability of ETFE/PTFE bottle adapters and ETFE/PTFE screw caps.

2.9 Storage conditions

Store the unit and accessories in a cool and dry place in cleaned condition only.

Storage temperature from -20 to 50°C (-4 to 122°F).

2.10 Dispenser Selection Chart

NOTICE Hydrofluoric acid (HF):

For dispensing hydrofluoric acid (HF), we recommend the Dispensette S Trace Analysis bottle-top dispenser with platinum-iridium valve spring. See separate operating instructions at www.brand.de/om.

	S	S Org.
1,2-diethylbenzene	✓	✓
1,4-dioxane	—	✓
1-butanol	✓	✓
1-decanol	✓	✓
Acetaldehyde	✓	✓
Acetone	✓	✓
Acetonitrile	✓	✓
Acetophenone	—	✓
Acetylacetone	✓	✓
Acetyl chloride	—	✓
Adipic acid	✓	—
Allyl alcohol	✓	✓
Aluminum chloride	✓	—
Formic acid, ≤ 100%	—	✓
Amino acids	✓	—
Ammonia, ≤ 20%	✓	✓
Ammonia, 20–30%	—	✓
Ammonium chloride	✓	—
Ammonium fluoride	✓	—
Ammonium sulfate	✓	—
Amyl alcohol (pentanol)	✓	✓
Amyl chloride (chloropentane)	—	✓
Aniline	✓	✓
Essential oils (aromatic oils)	—	✓
Barium chloride	✓	—
Benzaldehyde	✓	✓
Gasoline (petroleum spirit), boiling point 70–180°C	—	✓
Benzoic acid methyl ester	✓	✓
Benzene	✓	✓
Benzoyl chloride	✓	✓
Benzyl alcohol	✓	✓
Benzylamine	✓	✓
Benzyl chloride	✓	✓
Boric acid, ≤ 10%	✓	✓
Pyruvic acid	✓	✓
Bromobenzene	✓	✓
Bromonaphthalene	✓	✓
Hydrobromic acid	—	✓
Butanediol	✓	✓
Butyric acid	✓	✓
Butyl acetate	✓	✓
Butylamine	✓	✓
Butyl methyl ether	✓	✓
Calcium carbonate	✓	—
Calcium chloride	✓	—

	S	S Org.
Calcium hydroxide	✓	—
Calcium hypochlorite	✓	—
Chloroacetaldehyde, ≤ 45%	✓	✓
Chloroacetone	✓	✓
Chlorobenzene	✓	✓
Chlorobutane	✓	✓
Chloroacetic acid	✓	✓
Chloronaphtalene	✓	✓
Chloroform	—	✓
Chlorosulfonic acid	—	✓
Chromic acid, ≤ 50%	✓	✓
Chromosulfuric acid	✓	—
Cumene (isopropylbenzene)	✓	✓
Cyclohexane	—	✓
Cyclohexanone	✓	✓
Cyclopentane	—	✓
Decane	✓	✓
Dibenzyl ether	✓	✓
Dichlorobenzene	✓	✓
Dichloroacetic acid	—	✓
Dichloroethane	—	✓
Dichloroethylene	—	✓
Dichloromethane	—	✓
Diesel oil (fuel oil), boiling point 250–350°C	—	✓
Diethanolamine	✓	✓
Diethylamine	✓	✓
Diethylene glycol	✓	✓
Diethyl ether	—	✓
Dimethylaniline	✓	—
Dimethylformamide (DMF)	✓	✓
Dimethyl sulfoxide (DMSO)	✓	✓
Diphenyl ether	✓	✓
Acetic acid, ≤ 96%	✓	✓
Acetic acid, 100% (glacial acetic acid)	✓	✓
Acetic anhydride	—	✓
Ethanol	✓	✓
Ethanolamine	✓	✓
Ethyl acetate	✓	✓
Ethylbenzene	—	✓
Ethylene chloride	—	✓
Ethyl methyl ketone	✓	✓
Fluoroacetic acid	—	✓
Formaldehyde, ≤ 40%	✓	—
Formamide	✓	✓
Glycol (ethylene glycol)	✓	✓
Glycolic acid, ≤ 50%	✓	—
Glycerin	✓	✓
Urea	✓	—
Fuel oil (diesel oil), boiling point 250–350°C	—	✓
Heptane	—	✓
Hexane	—	✓
Hexanol	✓	✓
Hexanoic acid	✓	✓
Hydroiodic acid, ≤ 57% **	✓	✓

	S	S Org.
Isoamyl alcohol	✓	✓
Isobutanol	✓	✓
Isooctane	—	✓
Isopropanol (2-propanol)	✓	✓
Isopropyl ether	✓	✓
Potassium chloride	✓	—
Potassium dichromate	✓	—
Potassium hydroxide	✓	—
Potassium permanganate	✓	—
Cresol	—	✓
Copper sulfate	✓	—
Methanol	✓	✓
Methoxybenzene	✓	✓
Methyl tert-butyl ether (MTBE)	✓	✓
Methylene chloride	—	✓
Methyl formate	✓	✓
Methyl propyl ketone	✓	✓
Lactic acid	✓	—
Mineral oil (engine oil)	✓	✓
Monochloroacetic acid	✓	✓
n-amyl acetate	✓	✓
Sodium acetate	✓	—
Sodium chloride	✓	—
Sodium dichromate	✓	—
Sodium fluoride	✓	—
Sodium hypochlorite	✓	—
Sodium hydroxide, ≤ 30%	✓	—
Nitrobenzene	✓	✓
Oleic acid	✓	✓
Oxalic acid	✓	—
Pentane	—	✓
Perchloroethylene	—	✓
Perchloric acid	✓	✓
Peracetic acid	—	✓
Petroleum ether, boiling point 40–70°C	—	✓
Petroleum, boiling point 180–220°C	—	✓
Phenol	✓	✓
Phenylethanol	✓	✓
Phenyldiazine	✓	✓
Phosphoric acid, 85	✓	✓
Piperidine	✓	✓
Propionic acid	✓	✓
Propylene glycol (propanediol)	✓	✓
Pyridine	✓	✓
Salicylaldehyde	✓	✓
Nitric acid, ≤ 30%	✓	✓
Nitric acid, 30–70% */**/*	—	✓
Hydrochloric acid, ≤ 20%	✓	✓
Hydrochloric acid, 20–37% **	—	✓
Sulfuric acid, ≤ 98%	✓	✓
Silver acetate	✓	—
Silver nitrate	✓	—
Scintillation cocktail	✓	✓
Turpentine	—	✓

	S	S Org.
Tetrachloroethylene	—	✓
Carbon tetrachloride	—	✓
Tetrahydrofuran (THF) */**	—	✓
Tetramethylammonium hydroxide	✓	—
Toluene	—	✓
Trichlorobenzene	—	✓
Trichloroacetic acid	—	✓
Trichloroethane	—	✓
Trichloroethylene	—	✓
Trichlorotrifluoroethane	—	✓
Triethanolamine	✓	✓
Triethylene glycol	✓	✓
Trifluoroacetic acid (TFA)	—	✓
Trifluoroethane	—	✓
Hydrogen peroxide, ≤ 35%	—	✓
Tartaric acid	✓	—
Xylene	—	✓

	S	S Org.
Zinc chloride, ≤ 10%	✓	—
Zinc sulfate, ≤ 10%	✓	—

Legend:

✓ = The device is suitable for the medium

— = The device is not suitable for the medium

This table has been carefully tested and is based on the most current information available. Always observe the operating manual of the instrument and the specifications provided by the reagent manufacturer. In ad-

dition to the chemicals listed above, a number of organic or inorganic saline solutions (e.g. biological buffers), biological detergents as well as media for cell cultures can be dispensed. If you need chemical hazard statements that are not on the list, you are welcome to contact BRAND.

Last updated: 1024/16


* Use ETFE/PTFE bottle adapter

** Use PTFE sealing ring for valve block


*** Dispensette S Organic, Analog-adjustable.

3 Functional and operational components

Dispensette® S | Analog | with recirculation

Dispensette® S | Digital | with recirculation

4 Assembly

⚠ WARNING

Risk of injury from hazardous media

- Wear your personal protective equipment, in particular protective clothing suitable for the medium to be used, eye protection and protective gloves!
- Before use, find out how to use the medium safely and whether the laboratory instrument is suitable for this purpose.
- Observe the chemical-specific safety data sheets and follow the operating instructions for the medium at the workplace!

4.1 First Steps

1. Installing the filling tube/recirculation tube

- a. Adjust the length of the telescopic filling tube to the bottle height and attach it.

NOTICE Set the filling tube to approx. 2 cm above the bottom of the bottle so that no crystallized substances or solids are sucked into the instrument.

- b. Center the filling tube (side with smaller diameter) and attach it carefully to avoid damaging the olive-shaped nozzle.
- ⇒ If a dispensing cannula with recirculation valve is used, the recirculation tube must also be fitted.
- c. Insert the recirculation tube with the opening pointing outward.

2. Mounting and aligning the instrument on a bottle

- a. Screw the instrument (GL 45 threads) onto the reagent bottle and then align the dispensing cannula with the bottle label. This is done by rotating the valve block with the dispensing cannula.

NOTICE Do not tip

To prevent tipping, use a bottle stand – particularly with small bottles and when using the flexible discharge tube.

3. Transporting the instrument

⚠ WARNING! Hazardous media:

If hazardous media are used, handle the instrument and bottle only with protective gloves.

- a. For bottles with other thread sizes, select a suitable bottle adapter.
- b. When mounted to a reagent bottle, always carry the instrument as shown in the figure.

NOTICE Use the correct bottle adapters

The scope of delivery includes bottle adapters made of polypropylene (PP). These may be used only for media that do not attack PP. Alternatively, ETFE/PTFE bottle adapters can be used (Accessories/spare parts, p. 89). The user must review carefully to ensure the suitability of ETFE/PTFE bottle adapters.

NOTICE Template for adapter

An adapter template is available for download at www.brand.de in the Service section. This can be used to determine thread and adapter sizes. This adapter template for the respective product can also be found in the on-line shop <https://shop.brand.de> in the "Selection Charts" section.

4.2 Priming

NOTICE Rinse instrument thoroughly:

Before using the instrument for the first time, ensure it is rinsed carefully and discard the first few samples dispensed. Avoid splashes, as there may be ethanol and glycerine residues in the instrument.

Instruments with recirculation valve

- a. Open the screw cap of the discharge tube.

⚠ WARNING! Spraying of liquid:

Never press down the piston when the discharge tube is sealed with the screw cap!

Carefully open the screw cap, as liquid residues may collect there.

b. Turn the valve to 'Re-dose'.

c. For priming, pull the piston up approx. 30 mm and push it down until the lower stop. Repeat this procedure at least 5 times.

d. Turn the valve to 'Dispense'.

e. To avoid splashing, hold the opening of the dispensing cannula against the inside of a suitable collection container and dispense until the air has been removed from the dispensing cannula and no bubbles remain. Wipe away any remaining drops from the cannula.

Instruments without recirculation valve

- a.** Open the screw cap of the dispensing cannula (see 'Instruments with recirculation valve', Fig. a). To avoid splashes, hold dispensing cannula orifice on the inner wall of a suitable receiving vessel.
- b.** For priming pull the piston up approx. 30 mm and push it down until the lower stop. Repeat this process about 5 times until the air has been removed from the dispensing cannula and no bubbles remain.

5 Operation

⚠ WARNING


To be observed for every application, especially with hazardous media

- > Wear your personal protective equipment, in particular protective clothing, eye protection and protective gloves!
- > Before use, find out how to use the medium safely and whether the laboratory instrument is suitable for this purpose. Observe the safety data sheets and follow the operating instructions for the medium!
- > Never press down the piston when the dispensing cannula is sealed with the screw cap. Dispense slowly to avoid splashing.
- > Open the screw cap slowly to avoid splashing. Media residue can accumulate in the screw cap.

5.1 Dispensing

1. Select volume

Digital: Turn the volume adjustment knob until the desired volume is displayed (mechanical counter).

Analog: Loosen the volume selector thumb screw by a $\frac{3}{4}$ turn (1), slide the pointer vertically to the desired volume (2) and re-tighten the volume thumb screw (3).

Fix: The volume is permanently set and cannot be changed.

2. Dispensing

- a. Unscrew the screw cap of the dispensing cannula.
- b. For instruments with a recirculation valve, turn the valve to dispensing.
- c. Hold the dispensing cannula orifice on the inner wall of a suitable receiving vessel.

- d. Gently pull the piston up to the stop and then push it back down to the lower stop evenly and without applying too much force.

- e. Wipe off the dispensing cannula against the inner wall of the receiving vessel.
f. Seal the dispensing tube with the screw cap.

NOTICE Park position:

After use, always press the plunger down to the bottom stop (**park position**). If the piston is not pressed down to the bottom stop, media may escape unintentionally.

5.2 Accessories

5.2.1 Flexible discharge tube with recirculation valve

⚠ WARNING

Risk of media contact if the hose is damaged

If the hose is damaged, media contact may occur. Check the hose for damage (e.g. kinks and the like) before each use and replace it if necessary.


- Use a face shield and personal protective equipment.

- Secure the bottle with a bottle stand.
- In order to avoid splashing reagent, hold the dispensing tube firmly, and insert it into the holder provided after use.
- Rinse the tubing to clean it.
- Do not disassemble.

For serial dispensing, the flexible discharge tubing for the bottle-top dispensers Dispensette® S and Dispensette® S Organic can also be used (Accessories/spare parts, p. 89).


The specified accuracy and coefficient of variation of the instrument are only obtained for volumes > 2 ml and by gently approaching the upper and lower stops. The coil of the tubing can be stretched to a maximum length of 800 mm. Before use, make sure that the hose is properly looped and not twisted. The operating exclusions for the respective instrument in use apply.


Prerequisite:

- If the instrument was in use, it must be cleaned before installing the flexible discharge tube (Cleaning, p. 74).
- a. For instruments with a recirculation valve, set the valve to 'Recirculation' and pull the valve toggle upwards.
- b. Slide the discharge tube housing all the way up, then pull it forward with gentle up and down motions.

- c. Push the flexible discharge tube holder from the bottom of the valve block and tighten it. The Instrument must not be mounted on a bottle for this purpose. Fit the collecting tube.
- d. Press down the plug of the recirculation valve.


- e. Slide the flexible discharge tube housing into the valve block up to the stop.

- f. Push the housing all the way down.
- g. Put on the valve toggle made for the discharge valve and press it in firmly. In doing so, pay attention to the color coding and label (see assembly instructions 'Flexible discharge tubing for Dispensette® S').

NOTICE Use bottle stand, see accessories/spare parts, p. 89.

5.2.2 Drying tube

For moisture- or CO₂-sensitive media, it may be necessary to use a drying tube filled with suitable absorbent (not included in the scope of delivery).

See also Accessories/spare parts, p. 89.

- a. Unscrew the venting stopper using the installation tool.

- b. Screw in the filled drying tube.
- c. Place the PTFE sealing ring on the bottle thread or place the screwed-on bottle adapter and screw the Instrument onto the bottle.

NOTICE If necessary, seal the threads of the drying tube, the bottle and/or the bottle adapter with PTFE tape.

5.2.3 Sealing ring for valve block

For highly volatile or outgassing media, we recommend sealing the connection from the valve block to the bottle with the PTFE sealing ring and PTFE tape (Accessories/spare parts, p. 89).

5.2.4 Venting stopper for microfilter with Luer cone

For sterile media, we recommend the venting stopper with Luer cone for connecting a microfilter. This offers increased protection against contamination from the drawn-in air (Accessories/spare parts, p. 89).

- a. Unscrew the venting stopper (see 'Fitting the drying tube', Fig. a).
- b. Screw in the venting stopper with Luer cone.
- c. Place the PTFE sealing ring on the bottle thread or place the screwed-on bottle adapter and screw the Instrument onto the bottle.
- d. Insert a commercially available sterile filter into the Luer cone.

6 Error limits

Error limits in relation to the nominal capacity (= maximum volume) indicated on the instrument at equal temperature (20 °C/68 °F) of instrument, ambient environment and distilled water. The test was carried out in accordance with DIN EN ISO 8655-6 with a completely filled instrument and uniform and jerk-free dispensing.

Error limits

Nominal volume ml	A* ≤ ± %	µl	CV* ≤ %	µl
1	0.5	5	0.1	1
2	0.5	10	0.1	2
5	0.5	25	0.1	5
10	0.5	50	0.1	10
25	0.5	125	0.1	25
50	0.5	250	0.1	50
100	0.5	500	0.1	100

Partial volumes

The percentage values for A and CV are relative to the nominal volume (V_N) and must be converted for partial volumes (V_T).

$$A_T = \frac{V_N}{V_T} \cdot A_N$$

e.g.	Volume	A* ≤ ± %	µl	CV* ≤ %	µl
V_N	25.0	0.5	125	0.1	25
$V_T = 50\% N$	12.5	1.0	125	0.2	25
$V_T = 10\% N$	2.5	5.0	125	1.0	25
*A = Accuracy, CV = Coefficient of Variation					

NOTICE The error limits of DIN EN ISO 8655-5 are clearly undershot. The maximum total error for a single measurement can be calculated approximately from the sum of the error limits $FG = R + 2 VK$ (for the 25 ml size: $125 \mu l + 2 \times 25 \mu l = 175 \mu l$).

7 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3–12 months. This time frame should be adjusted to correspond with individual requirements.

Gravimetric volume testing according to DIN EN ISO 8655- 6 (for measurement conditions, see Error limits, p. 69) is performed as follows:

1. Prepare instrument

Clean the instrument (Cleaning, p. 74), fill it with test liquid (distilled water) and then prime it carefully.

2. Check volume

- Carry out 10 doses with test liquid in 3 volume ranges (100 %, 50 %, 10 %).
- For discharge, depress the piston slowly and steadily without force until the lower stop
- Wipe off the tip of the dispensing tube.
- Weigh the dispensed quantity on an analytical balance. Refer to the operating manual of the balance manufacturer.
- Calculate the dispensed volume. The Z factor takes account of the temperature and air buoyancy.

Calculation (for nominal volume)

x_i = weighing results

n = number of weighings

V_0 = nominal volume

Z = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPa)

Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \bar{x} * Z$$

Accuracy*:

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Standard deviation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

Coefficient of variation*:

$$CV\% = \frac{100 s}{\bar{V}}$$

*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

NOTICE Test Instructions (SOPs) are available for download from www.brand.de.

8 Adjustment

After a long period of usage, an adjustment of the instrument might be necessary.

- Calibrate for example at nominal volume (Checking the Volume (Calibration)).
- Calculate mean volume (actual value) (Checking the Volume (Calibration)).
- Adjust the instrument (to the calculated mean volume).
- After the adjustment, further calibration is necessary to confirm the appropriate adjustment.

Example

The gravimetric control of a 10 ml instrument results in a current value of 9.90 ml with a set nominal volume of 10.00 ml.

8.1 Digital model


- Push the lock to the left and remove the front casing.


- Pull out the safety slider. The calibration cover will come off in the process. Dispose of the calibration cover.


- Pull out the red rotary knob to decouple the gears. Set the determined actual value (e.g. 9.90 ml).

Nominal volume [ml]	Digital max. +/- [μ l]
1	12
2	24
5	60
10	120
25	300
50	600

- d.** Push the red rotary knob and then the safety slider back in.


- e.** Close the casing and push the lock to the right.
- ⇒ The calibration is complete.
 - ⇒ The change in factory calibration is indicated by a red marking (circle).

8.2 Analog-adjustable model

- a.** Insert the pin of the installation tool into the calibration cover (pos. 1) and break it off by turning it. Dispose of the calibration cover.


- b.** Insert the pin of the installation tool into the calibration screw and turn to the left to increase the dispensing volume or turn to the right to decrease the dispensing volume (e.g. current value 9.97 ml approx. 1/2 turn to the left).

Nominal volume	Analog/Fixed max. +/- [μl]	One rotation corresponds to [μl]
1 ml	6	~ 8
2 ml	12	~ 16
5 ml	30	~ 40
10 ml	60	~ 80
25 ml	150	~ 130
50 ml	300	~ 265
100 ml	600	~ 400

c. The calibration is complete.

⇒ The change in calibration is indicated by a red disk (circle in illustration).

9 Cleaning

⚠ WARNING

There is a risk of exposure to liquids during cleaning.

The components cylinder, valves, suction pipe and dosing cannula are filled with liquid!

- Wear protective clothing, eye protection and protective gloves to protect yourself from hazards caused by the liquid during cleaning.
- Never point the openings of the suction tube, dosing cannula and valves towards the body in order to avoid hazards from escaping liquid.
- Never remove the dosing cannula when the dosing cylinder is full, as liquid will then escape.

The instrument must be cleaned in the following situations to assure correct operation:

- Before using the instrument for the first time.
- Immediately if the piston becomes stiff.
- Before you change the liquid.
- Before storing the instrument.
- After you have stored the instrument.
- Before dismantling the instrument.
- Before autoclaving the instrument.
- Before replacing one or more valves.
- Regularly when using liquids which form deposits (e.g., crystallizing liquids).
- Regularly if liquids accumulate in the screw cap
- Regularly if you use hygroscopic liquids (e.g. concentrated sulphuric acid).

9.1 Cleaning type Analog, type Fixed

1. Empty the instrument completely

- a. Screw the instrument onto an empty bottle and empty it completely by dispensing. If the instrument is equipped with a recirculation valve, it must be emptied in the dosing and recirculation position.

2. Rinse instrument

- a. Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and fill and empty completely several times to rinse.

3. Dismantling the piston


- a. Hold the housing shells firmly and unscrew the piston mount completely by turning it to the left.

NOTICE Do not interchange pistons

The pistons of the instruments are individually fitted and must not be interchanged with pistons from other instruments!

- b. Carefully pull out the piston.

4. Clean the piston and cylinder, and refit the instrument

- a. Clean the piston and cylinder. Carefully remove any deposits from the upper edge of the dosing cylinder (e.g., with a bottle brush and thin soapy water).
- b. Rinse the piston and cylinder with deionized water, and dry carefully.
- c. Attach the dosing cannula to open the discharge valve.

d.

- d. Grasp the piston at the tip and carefully push it fully into the cylinder **vertically while rotating it** and re-assemble the instrument.

9.2 Cleaning Digital model

1. Empty the instrument completely

- a. Screw the instrument onto an empty bottle and empty it completely by dispensing. If the instrument is equipped with a recirculation valve, it must be emptied in the dosing and recirculation position.

2. Rinse instrument

- a. Screw the instrument onto a bottle filled with a suitable cleaning agent (e.g. deionized water) and fill and empty completely several times to rinse.

3. Disassemble piston

a.


- a. Assembly and disassembly must only be carried out at the set maximum volume. To do this, turn the volume adjustment knob in the + direction to the stop (= maximum volume).

NOTICE Do not interchange pistons

The pistons of the instruments are individually fitted and must not be interchanged with pistons from other instruments!


- b. Push the lock to the left.
- c. Remove the front casing.

- d. Loosen the piston mounting nut with the installation tool and pull the rear casing part out to the back.


- e. Then carefully pull the piston out of the cylinder.

4. Clean the piston and cylinder, and refit the instrument

- a. Clean the piston and the cylinder. Carefully remove any deposits from the upper edge of the dispensing cylinder (e.g., with a bottle brush and thin soapy water).
- b. Rinse the piston and cylinder with deionized water, and dry carefully.
- c. Attach the discharge tube to open the discharge valve.


- d. Grasp the piston at the tip and carefully push it fully into the cylinder **vertically while rotating it** and re-assemble the instrument.


NOTICE Dispensette S/S Organic Digital:

Correctly position the stop segment during assembly.

The stop segment must engage under the stop ring of the cylinder. When tightening the piston mounting nut with the installation tool, the piston/cylinder unit must be pressed backward toward the rear casing using your thumb.

Correct


Incorrect

9.3 Replacing the dosing cannula/valves


NOTICE A functional test must always be carried out after replacing components.

9.3.1 Replacing the dosing cannula


- a. For instruments with a recirculation valve, set the valve to 'Recirculation'.


- b. Pull the valve toggle upwards.

- c. Slide the dispensing cannula housing all the way up, then pull it forward with gentle up and down motions.

- d.** Slightly push the dosing cannula onto the valve.

Then pull the dosing cannula slightly upwards.
The coupling piece slides downwards

Push the dosing cannula onto the valve as far as it will go.
Press the dosing cannula down to lock it.

- e.** Push the housing all the the discharge tube way down.

- f.** For units with a redosing valve, place the valve knob in the 'redosing' position and push it downwards.

- g.** Carry out a functional test after replacement.

9.3.2 Replacing the valves

NOTICE Use the correct valves

Always install the valve intended for the particular instrument model and size! (See ordering data, Ordering Information, p. 87).

Identical filling valves but different discharge valves are used for Dispensette® S and Dispensette® S Organic. For easy identification, the discharge valves for the Dispensette® S Organic are marked with 'ORG'.

9.3.2.1 Replacing the discharge valve

- a. After disassembling the dispensing cannula (see Replacing the dosing cannula, p. 79), use the assembly tool to unscrew the discharge valve.

- b. Firmly screw in the new discharge valve first by hand and then tighten it with the assembly tool. The thread should not be visible.

NOTICE Safety ball in the discharge valve


The discharge valve has a safety ball that closes automatically when the discharge tube is removed. Attaching the discharge tube reopens the safety ball.

9.3.2.2 Replacing the filling valve

- a. Pull out the recirculation tube and the telescopic filling tube.

- b. Use the assembly tool to unscrew the filling valve.
- c. Screw in the new filling valve first by hand and then tighten it with the assembly tool.

9.3.2.3 Loosen stuck filling valve

If the instrument does not fill up, and if some elastic resistance is evident when the piston is pulled upward, then it is possible that the valve ball is stuck.

In this case, dislodge the valve ball using light pressure; for example, with a 200 µl plastic pipette tip (see the figure at the side).

10 Autoclaving

The instrument can be autoclaved at 121°C (250°F), 2 bar, and a holding time of at least 15 minutes according to DIN EN 285.

- a. The instrument must be carefully cleaned prior to autoclaving (Cleaning, p. 74).
- b. Open the screw cap on the discharge tube, and for instruments with a recirculation valve, set the valve to "Dispense".
- c. For all instruments:
Check the tight fit of the filling valve (Fig. 1a) and, if necessary, tighten the filling valve prior to autoclaving.

For the Digital type:

Additionally check the tight fit of the piston mounting (Fig. 1b) and, if necessary, tighten the piston mounting before autoclaving.



- d. To ensure unhindered access for the steam and to prevent the valve ball in the filling valve from possibly becoming stuck, hold the instrument with the discharge piston pressed vertically downward, and gently tap against the casing with your hand (Fig. 2).

Then lay it horizontally in the autoclave.

NOTICE Be sure to avoid the instrument coming into contact with metal surfaces in the autoclave!

NOTICE Allow the instrument to cool down

Do not reassemble the instrument until it has cooled down to room temperature (Cooling time approx. 2 hours). After every autoclaving, inspect all parts for deformities or damage. If necessary, replace them. The effectiveness of autoclaving must be verified by the user.

11 Troubleshooting

Problem	Possible cause	Corrective action
Piston sluggish or stuck.	Crystal deposits Contaminants	Stop dispensing immediately. Loosen piston with circular motion, but do not disassemble. Perform cleaning (Cleaning, p. 74).
	The dispensing cylinder or the piston is broken.	If the dispensing cylinder or piston is broken, replace the instrument. A repair is not possible in this case.
Filling not possible.	Volume adjusted to minimum setting.	Set the desired volume (Dispensing, p. 65).
	Filling valve stuck.	Unscrew the filling valve from the valve block and clean it. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it (Loosen stuck filling valve, p. 82). If necessary, replace the filling valve.
Dispensing not possible.	Discharge valve stuck.	Unscrew the discharge valve from the valve block, clean it and replace if necessary. If the valve ball is stuck, use a 200 µl plastic tip to dislodge it.
Discharge tube (with or without recirculation valve) cannot be mounted sufficiently.	Discharge valve is not screwed in deeply enough.	Tighten the discharge valve with the installation tool until it meets the stop so that the threads are no longer visible.
Air bubbles in the instrument.	Reagent with high vapor pressure has been drawn in too quickly.	Slowly draw in reagent.
	Valve screw connections loose.	Tighten the valves firmly with the installation tool.
	The instrument has not been primed.	Prime the instrument (Priming, p. 63).
	Filling tube is loose or damaged.	Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube.
	Valves not firmly connected or damaged.	Perform cleaning (Cleaning, p. 74). Tighten the valves using the installation tool.
Dispensed volume too low.	Filling tube is loose or damaged.	Perform cleaning (Cleaning, p. 74). Push the filling tube on firmly. If necessary cut off approx. 1 cm of tube at the upper end and re-connect it or replace filling tube (Replacing the dosing cannula, p. 79).
	Filling valve is loose or damaged.	Perform cleaning (Cleaning, p. 74). Tighten the filling valve using the installation tool. If necessary, replace filling valve.

Problem	Possible cause	Corrective action
Leaking liquid between instrument and bottle.	Recirculation tube not connected.	Connect recirculation tube (First Steps, p. 62, Fig. 3).
	Volatile reagent dispensed without sealing ring for valve block.	Install sealing ring for valve block (First Steps, p. 62).
Liquid discharge at the upper edge of the dispensing cylinder.	Piston not in lowest position after dispensing.	After use, always depress the piston to the lower stop (put it in the “parking position”) (see Dispensing, p. 65)
	Hygroscopic dispensing medium leads to a leakage of medium.	Regularly wipe off any accumulation of liquid on the upper edge of the dispensing cylinder.

12 Product markings

Symbol or number	Meaning
	General warning symbol
	Follow the operating instructions
	Use eye protection
	Use hand protection
	Wear protective clothing
XXZXXXXX	Serial number
 25	The instrument is marked in accordance with the German Measurement and Calibration Act as well as the Measurement and Calibration Regulation. DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was affixed.
www.brand.de/ip	Patent information

13 Ordering Information

Dispensette® S, Digital

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.1 – 1	0.005	4600310	4600311
0.2 – 2	0.01	4600320	4600321
0.5 – 5	0.02	4600330	4600331
1 – 10	0.05	4600340	4600341
2.5 - 25	0.1	4600350	4600351
5 – 50	0.2	4600360	4600361

Dispensette® S, Analog-adjustable

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.1 - 1	0.02	4600100	4600101
0.2 - 2	0.05	4600120	4600121
0.5 - 5	0.1	4600130	4600131
1 - 10	0.2	4600140	4600141
2.5 - 25	0.5	4600150	4600151
5 - 50	1.0	4600160	4600161
10 - 100	1.0	4600170	4600171

Dispensette® S, Fixed-volume

Volume ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
1	4600210	4600211
2	4600220	4600221
5	4600230	4600231
10	4600240	4600241
Fixed volumes by request: 0.5–100 ml (please specify when ordering)	4600290	4600291

NOTICE see Scope of delivery, p. 52

Dispensette® S Organic, Digital

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.5 - 5	0.02	4630330	4630331
1 - 10	0.05	4630340	4630341
2.5 - 25	0.1	4630350	4630351
5 - 50	0.2	4630360	4630361

Dispensette® S Organic, Analog-adjustable

Volume ml	Graduation ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
0.5 - 5	0.1	4630130	4630131
1 - 10	0.2	4630140	4630141
2.5 - 25	0.5	4630150	4630151
5 - 50	1.0	4630160	4630161
10 - 100	1.0	4630170	4630171

Dispensette® S Organic, Fixed-volume

Volume ml	Discharge tube without recirculation valve Order no.	Discharge tube with recirculation valve Order no.
5	4630230	4630231
10	4630240	4630241
Fixed volumes by request: 2–100 ml (please specify when ordering)	4630290	4630291

NOTICE Dispensing hydrofluoric acid (HF): For dispensing hydrofluoric acid (HF), we recommend the bottle-top dispenser Dispensette® S Trace Analysis with platinum-iridium valve spring (see separate operating instructions).

14 Accessories/spare parts

Bottle adapter

External thread	for bottle thread/ground joint	Material	Order No.
GL 32	GL 25	PP	704325
GL 32	GL 28/ S 28	PP	704328
GL 32	GL 30	PP	704330
GL 32	GL 45	PP	704345
GL 45	GL 32	PP	704396
GL 45	GL 35	PP	704431
GL 45	GL 38	PP	704397
GL 45	S* 40	PP	704343
GL 45	S* 50	PP	704350
GL 45	S* 54	PP	704430
GL 45	S* 60	PP	704348
GL 32	GL 25	ETFE	704375
GL 32	GL 28/ S 28	ETFE	704378
GL 32	GL 45	ETFE	704395
GL 45	GL 32	ETFE	704398
GL 45	GL 38	ETFE	704399
GL 45	S* 40	PTFE	704391
GL 32	NS 19/26	PP	704419
GL 32	NS 24/29	PP	704424
GL 32	NS 29/32	PP	704429

* buttress thread

NOTICE Template for adapter:

An adapter template is available for download at www.brand.de in the Service section. This can be used to determine thread and adapter sizes. This adapter template for the respective product can also be found in the on-line shop store.brand.de in the "Selection Charts" section.

Discharge tubes without recirculation valve

Discharge tube without recirculation valve for Dispensette® S

Packaging unit Pack of 1.

Nominal volume ml	Model	Length (mm)	Order no.
1, 2, 5, 10	Fine tip	108	708002
5, 10	Standard	108	708005
25, 50, 100	Fine tip	135	708006
25, 50, 100	Standard	135	708008

Discharge tubes without recirculation valve for Dispensette® S Organic

Packaging unit Pack of 1.

Nominal volume ml	Model	Length (mm)	Order no.
5, 10	Fine tip	108	708012
5, 10	Standard	108	708014
25, 50, 100	Fine tip	135	708016
25, 50, 100	Standard	135	708019

Discharge tubes with recirculation valve

Discharge tubes with recirculation valve for Dispensette® S

Packaging unit Pack of 1.

Nominal volume ml	Model	Length (mm)	Order no.
1, 2, 5, 10	Fine tip	108	708102
5, 10	Standard	108	708104
25, 50, 100	Fine tip	135	708106
25, 50, 100	Standard	135	708109

Discharge tubes with recirculation valve for Dispensette® S Organic

Packaging unit Pack of 1.

Nominal volume ml	Model	Length (mm)	Order no.
5, 10	Fine tip	108	708112
5, 10	Standard	108	708114
25, 50, 100	Fine tip	135	708116
25, 50, 100	Standard	135	708119

Flexible discharge tubing

PTFE, spiraled, approx. 800 mm long with safety handle.

Packaging unit Pack of 1.

Not suitable for hydrofluoric acid (HF)

Nominal volume ml	Discharge tubing outer diameter	Discharge tubing inner diameter	Order no.
1, 2, 5, 10	3	2	708132
25, 50, 100	4.5	3	708134

Discharge valve Dispensette® S

PFA/Boro 3.3/ceramic/platinum-iridium.

No valve marking.

Pkg unit Pack of 1.

for nominal volume ml	Order no.
1, 2*	6749
5, 10	6727
25, 50, 100	6728

* with valve marking '1 + 2'

Dispensette® S Organic discharge valve

PFA/Boro 3.3/ceramic/tantalum.

Valve marked with 'ORG'.

Pkg unit Pack of 1.

for nominal volume ml	Order no.
5, 10	6729
25, 50, 100	6730

Dispensette® S and Dispensette® S Organic filling valve

PFA/ETFE/Boro 3.3/ceramic.


No valve marking.

Pkg unit Pack of 1.

for nominal volume ml	Order no.
1, 2, 5, 10	6734
25, 50, 100	6735

Telescopic filling tubes for Dispensette® S and Dispensette® S Organic

FEP. Individually adjustable lengths.

Packaging unit Pack of 1.


for nominal volume ml	Outer diameter mm	Length (mm)	Order no.
1, 2, 5, 10	6	70-140	708210
		125-240	708212
		195-350	708214
		250-480	708216
25, 50, 100	7.6	170-330	708218
		250-480	708220

Venting stopper for microfilter with Luer cone

Description	Packaging unit	Order no.
Venting stopper for microfilter with Luer cone. PP. Venting stopper and PTFE-sealing ring.	1 pc.	704495

Recirculation tube

Description	Packaging unit	Order no.
Recirculation tube. FEP	1 pc.	6747

Adjustment/installation tool

Description	Packaging unit	Order no.
Adjustment/installation tool	1 pc.	6748

Sealing ring for valve block

Description	Packaging unit	Order no.
Sealing ring for valve block. PTFE, for highly volatile media.	1 pc.	704486

Bottle stand

Description	Packaging unit	Order no.
Bottle stand PP. Support rod, 325 mm, base plate 220 x 160 mm.	1 pc.	704275

Drying tube

Description	Packaging unit	Order no.
Drying tube, including PTFE sealing ring (without drying agent)	1 pc.	707930


Screw cap with fastener

Packaging unit Pack of 1.

Description	Nominal volume ml	Order No.
PP, red, for Dispensette® S	1, 2, 5, 10	706018
	25, 50, 100	706019
PP, yellow, for Dispensette® S Organic	5, 10	706025
	25, 50, 100	706027
ETFE, for Dispensette® S and Dispensette® S Organic	1, 2, 5, 10	706029
PTFE, for Dispensette® S and Dispensette® S Organic	25, 50, 100	706031

Dispensette® drum extraction system

Pkg unit Pack of 1.

Description	Order no.
Dispensette® drum extraction system	704261
Tripod connection for wall mount, stainless steel, Dispensette® extraction system	704268
Table/rack clamp for wall mount for Dispensette® extraction system	704272

15 Repairs

15.1 Sending for repair

NOTICE Transporting hazardous materials without approval is prohibited by law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in
www.brand.co.in

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562
F +1 - 860 - 767 2563
info@brandtech.com
www.brandtech.com

China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn
www.brand.cn.com

16 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the device. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive Liquids.

The complete SOP for testing can be downloaded from www.brand.de or www.brandtech.com.

BRAND also offers you the option of having your devices calibrated through our factory calibration service or through our accredited calibration laboratory. Just send us the devices to be calibrated, indicating the type of calibration you would like. You will get your devices back in a few days. A detailed calibration report (factory calibration) or an accredited calibration certificate in accordance with DIN EN ISO/IEC 17025 is enclosed with each device. More information can be obtained from your retailer or directly from BRAND. The order document is available for download at www.brand.de (Service & Support).

For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. Our service partners can forward your devices to BRAND for factory calibration, if required.

17 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

18 Disposal

Before disposal, observe the relevant national disposal regulations and ensure that the product is disposed of properly.

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