Betriebsanleitung
Operating instructions
Notice d'instructions
Istruzioni per l'uso

Hei-MIX Titramax 100/101/1000
Hei-MIX Vibramax 100
Hei-MIX Rotamax 120

 **heidolph**
research made easy

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Seite 4 – 23

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Contents

Introduction

About this document.....	28
Typographic conventions.....	28
Copyright protection.....	28
Safety instructions in the official languages of the European Union	28

General product information

Directives applied, product certification.....	29
Copyright notice	29
California Residents.....	29
Residual risk.....	29
Intended use.....	29
Compliant use.....	30
Reasonably foreseeable misuse	30
Transportation.....	30
Storage.....	30
Acclimatization	30
Permissible ambient conditions	30
Energy saving potential and cost efficiency	31

Safety

General safety instructions.....	31
Electrical safety.....	31
Operational safety	31
Occupational safety	32
Electrostatic charging of liquid or powdered media.....	32
Personal protective equipment (PPE).....	32
Environmental protection.....	32
Biohazard.....	32
Other regulations.....	33

Device description

Mechanical design, control elements	34
Movement type / shaking orbit	34
Commissioning	35
Positioning the device	35
Power supply	35
Placing and attaching laboratory vessels	36
Hei-MIX Titramax 100 / 101 / 1000	36
Hei-mix Vibramax 100 / Rotamax 120 / attachment with clamping rollers	37
Operation	38
Switching the device on/off	38
Operating modes	38

Troubleshooting

Troubleshooting	39
-----------------------	----

Appendix

Technical specifications	40
Performance range/maximum load	42
Scope of delivery	43
Accessories	43
Device service	44
General cleaning instructions	44
Replace the device fuse	44
Repairs - Return of equipment	45
Maintenance	45
Disposal	45
Warranty statement	46
Contact details	46
Certificate of decontamination	47

About this document

This operating manual describes all functions and operation of the Hei-MIX Titramax 100/101/1000, Hei-MIX Vibramax 100 and Hei-MIX Rotamax 120 platform shakers.

The operating instructions are an integral part of the described device!

Typographic conventions

This document uses the following symbols, signal words, and highlights:

	Description
	<p>Warning symbols in combination with a signal word indicate dangers:</p> <p>DANGER Indicates an immediate dangerous situation. Failure to respect the indications will result in death or serious injury.</p> <p>WARNING Indicates a potential danger. Failure to respect the indications will result in serious injuries.</p> <p>CAUTION Indicates a potential hazard. If not avoided, damage to property and minor to moderate injuries can occur.</p>
	<p>Mandatory signs are used to indicate important and useful information on handling a product. This information is used to ensure operational safety and to maintain the value of the product.</p>
[...]	<p>Square brackets identify labels of control elements on the device as well as labels and entries in software masks and user interfaces.</p>
→	<p>The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.</p>

Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, as well as utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Scientific Products GmbH. Any violation is subject to compensation for damage.

Safety instructions in the official languages of the European Union

A summary of all safety instructions in the official languages of the European Union can be found in our Safety Guide for the product group Hei-SHAKE (shakers and mixers, Ref. 01-005-006-99). This document is available for download on our homepage in the most up-to-date version.

Directives applied, product certification

CE marking

The devices comply with the following standards:

- 2014/35/EU, Low Voltage Directive
- EMC Directive 2004/108/EC
- RoHS 2 directive 2011/65/EU, including amendment 2015/863

Copyright notice

The software implemented in this product is protected by copyright laws. The rights holder is Heidolph Scientific Products GmbH, Walpersdorfer Straße 12, 91126 Schwabach/Germany. Any open source software components in this software are excluded from our copyright. Further information is available in the service area on our website www.heidolph.com.

Violations of copyright (e.g. unauthorized use or modification of the software) may result in civil claims (e.g. omission, compensation for damages) and/or entail criminal penalties within the meaning of the legislation of the destination country.

California Residents

Important information for California residents regarding Prop 65. Please visit www.P65Warnings.ca.gov for more information.

Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

Intended use

The described platform shaker in this document has been specially developed for the following tasks:

- Shaking
- Mixing
- Emulsifying
- Suspending
- Separating
- Dissolving
- Staining

The areas of application of the described platform shaker include chemical, biological and environmental-analytical laboratory and research applications, applications of basic research and comparable facilities.

Due to its design, the device in its delivery condition may only be used in analytical processes or in laboratory-like conditions in the food, cosmetics, and pharmaceutical industries as well as other comparable industries that manufacture products intended for consumption by humans or animals, or for use on humans or animals.

Any other use of this device is not considered as intended!

Compliant use

The compliance of each individual application must always be evaluated by the user. Any additional measures necessary to ensure compliance are always the responsibility of the user.

Reasonably foreseeable misuse

Additional measures may be necessary for use under conditions or for purposes deviating from the intended use, and/or specific guidelines and safety regulations must be observed. Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is the sole responsibility of the operator.

All risks resulting from improper use are borne solely by the operator.

The device may only be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the sole responsibility of the operator!

Transportation

During transportation, avoid severe shocks and mechanical stresses that can cause damage to the device. Keep the original packaging in a dry and protected place for later use.

Storage

Always store the product in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

Acclimatization

After each transportation and after storage under critical climatic conditions (e.g., high temperature difference between inside and outside), allow the device to acclimatize at room temperature for at least two hours to prevent possible damage due to condensation before commissioning it in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections only after the device has been acclimatized!

Permissible ambient conditions

The device is designed for indoor use only. The device is **NOT** suitable for outdoor use! The device is **NOT** suitable for use in potentially explosive areas!

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

Energy saving potential and cost efficiency

Even in standby mode, small amounts of energy are used to supply individual device functions. However, especially over longer periods of nonuse, this resource consumption adds up and causes avoidable costs.

Therefore, always switch off the device when not in use, especially for long periods of time (devices with main switch: Main switch in position [O]; devices without main switch: Unplug) to effectively avoid unnecessary energy consumption and the resulting costs.

General safety instructions

Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.

Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself and, where applicable, on connected devices or the supply connections.

If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.

Only use the device in accordance with the regulations for intended use.

Electrical safety

Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.

Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).

Always use the power supply cord provided with the device.

Prior to use, check that the device and the power supply cord are free of visible damage.

Make sure that the power plug can be reached directly at any time, in order to disconnect the device from the power supply at any time without delay.

Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Scientific Products GmbH.

Always disconnect the device from the power supply system before carrying out any maintenance, cleaning or repair work.

Operational safety

Operate the device under a closed ventilated fume hood when working with potentially hazardous substances (see EN 14175 and DIN 12924).

Do not make any unauthorized changes or modifications to the device!

Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!

Rectify malfunctions or faults on the device immediately.

Switch off and disconnect the device from the power supply, preventing reconnection, if it is not possible to eliminate the malfunction or rectify the fault immediately.

Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

Occupational safety

Always use the prescribed personal protective equipment (PPE) such as protective clothing, safety goggles, protective gloves, safety shoes, etc.

Do not operate any other devices in the immediate vicinity of the device ...

- which can generate electromagnetic fields in the frequency range between 9×10^3 Hz to 3×10^{11} Hz,
- which generate emission or radiation sources in the frequency range 3×10^{11} Hz to 3×10^{15} Hz (in the optical spectral range wavelengths from 1,000 μm to 0.1 μm),
- which generate ultrasonic or ionizing waves.

Do not process any substances that could release energy in an uncontrolled manner (e.g. self-ignition).

Do not process substances in which the energy input through mixing poses a danger.

Wipe off any liquid that may have spilled on the device immediately.

Electrostatic charging of liquid or powdered media

Liquid or powdered media as well as the laboratory vessels used can become dangerously electrostatically charged due to the shaking movement or when emptying or cleaning the vessel. The intensity of the charge generated depends on the properties of the medium, the flow rate, the process parameters as well as the size, geometry and material of the laboratory vessels (source: TRGS 727, Section 4).

The evaluation of a potential explosion hazard due to electrostatic charging and the definition of suitable safety measures are always the responsibility of the operator and are always based on the substances to be processed and the individual application.

Personal protective equipment (PPE)

The operating company must determine and provide the necessary PPE, depending on the respective application and the media and chemicals used.

The required instruction of the personnel is solely within the operating company's responsibility.

Environmental protection

When processing environmentally hazardous substances, take appropriate measures to avoid hazards to the environment.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking of the device with the biohazard warning symbol.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Scientific Products GmbH.

The operating company is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Mechanical design, control elements

The following illustration shows the mechanical structure of a Hei-MIX Titramax 100 platform shaker. The Hei-MIX Titramax101/1000, Vibramax 100 and Rotamax 120 units have identical controls.

Designation

- 1 Shaking platform
- 2 Rotary control [Timer]
- 3 Rotary control [Shaking intensity]
- 4 Selector switch [Timer]
- 5 On/off switch

Function

Surface for placing laboratory vessels; different variants are available, see scope of delivery and accessories.

Rotary control for preselecting the duration for continuous shaking (up to a maximum of 120 minutes).

Rotary control for preselecting the shaking intensity (setting range depends on the device, see technical data).

Selector for switching between continuous and timer operation:

- Position [0]: Unit operates in continuous operation
- Position [Timer]: Unit operates in timer mode

Main switch for switching the device on and off. When the device is switched on, the on/off button lights up green.

Movement type / shaking orbit

Model	Motion	Shaking Orbit
Vibramax 100	circular vibrating	3 mm
Titramax 100		1.5 mm
Titramax 101		3 mm
Titramax 1000		1.5 mm
Rotamax 120	orbital	20 mm

Positioning the device

CAUTION

Improper positioning of the device may result in direct and indirect damage to property due to falling and/or spilling of fluids!

The surface must be even, non-slip, proper and sufficiently stable.

The device must be freely accessible at all times.

Make sure that a sufficient safety distance to moving components can be maintained on all sides.

Observe the generally valid rules for the safe positioning of work equipment.

The professional mounting and positioning of the device including the supplied accessories is within the sole responsibility of the operator!

Heidolph Scientific Products GmbH shall under no circumstances be liable for any personal injury or damage to property resulting from failure to observe the instructions for the correct installation and positioning of the device.

Power supply

DANGER

Observe the instructions in section "Allgemeine Sicherheitshinweise" on page 7, in particular Electrical safety.

Connect the power supply cord

Before connecting the power supply cord, make sure that the main switch of the device is in position [0] (off).

Connect the cable coupling of the coded power supply cord to the IEC appliance inlet on the back of the device.

Connect the power supply cord to a properly secured mains socket-outlet.

Disconnecting the power supply cord

Before disconnecting the power supply cord, switch off the device.

Disconnect the power supply cord from the power socket.

Disconnect the power supply cord from the IEC appliance inlet on the back of the device.

Placing and attaching laboratory vessels

For placing and attaching your laboratory vessels, use the various Heidolph attachments which are available as an option. For detailed information on the available accessories, see www.heidolph.com → Products → Platform shakers → Accessories, as well as the Heidolph Assembly Instructions Hei-SHAKE, Ref. 01-005-007-13.

- At low shaking frequencies, vessels such as Petri dishes can also be placed without using any kind of attachment on the non-slip contoured rubber mat.
- Always place individual vessels in the center of the shaking platform.
 - Always distribute several vessels evenly on the shaking platform.

Hei-MIX Titramax 100 / 101 / 1000

Insert multi-well plates into the pockets of the contoured rubber mat supplied. The multi-well plates are securely fixed in the pockets.

Hei-mix Vibramax 100 / Rotamax 120 / attachment with clamping rollers

The attachment with clamping rollers (accessory) is designed to fix laboratory vessels on the shaking platform of the device.

Note that only vessels of the same diameter can be fixed in parallel!

Mount the attachment to the shaking platform, using the thumbscrews supplied.

Place the vessels on the non-slip contoured rubber mat.

Fix the vessels in parallel with the clamping rollers.

Switching the device on/off

To switch on and off, use the main switch on the front of of the device.

Operating modes

The platform shakers described can be operated in continuous or timer mode. Regardless of the operating mode selected, observe the following specific safety instructions.

WARNING

- Whenever possible, use closed vessels for processing corrosive, toxic or biohazardous substances and seal them safely.
- Increase the speed only gradually to the desired shaking intensity, especially with open vessels, and observe the fluid movements.
- Please note the maximum lifting capacity of the platform of your device, see section "" on page 40.
- Lower the shaking intensity or reduce the total load on the platform when the device starts to vibrate during operation.
- Always ensure that the vessels are evenly distributed on the platform.
- Pay particular attention to an even distribution of the weight load when vessels of different sizes and/or differently filled vessels are placed on the platform at the same time!
- Before switching on the device, make sure that all vessels are properly fixed on the platform.
- Use suitable attachments for the used vessels. Further information on available accessories can be found on our website at www.heidolph.com.
- Depending on the effort and the duration of use, there is a risk of over-exposure to vibration of your hand and arm when using the hand to hold the vessel during operation. Observe the European Directive 2002/44/EC regarding the permissible hand-arm vibration exposure.
- Always use the necessary and appropriate personal protective equipment!

Continuous operation

In continuous operation, the shaking movement starts as soon as the device is switched on via the main switch.

- Place the vessel(s) on the shaking platform.
- Switch on the device.
- Adjust the shaking frequency with the rotary knob 'shaking frequency adjustment'.
- The shaking movement stops as soon as the device is switched off.

Timer mode

In timer operation, a specific period of time (up to 120 minutes) can be set for continuous shaking.

- Place the vessel(s) on the shaking platform.
- Use the 'time adjustment' rotary knob to set a desired time for shaking operation.
- Switch on the device.
- Adjust the shaking frequency with the rotary knob 'shaking frequency adjustment'.
- The shaking movement stops as soon as the specified time has elapsed.

The shaking movement can be stopped manually at any time in timer mode if required. To do this, switch off the device using the main switch.

Troubleshooting

Failure	Cause	Possible remedy
Display remains dark after switching on	No mains voltage	Check the power supply cord for damages/correct connection
	Power supply cord defective	Check building fuses
LED of the [Start/Stop] button does not light up when the function is activated	Display defective	Contact Heidolph Sales or Service
	Button LED defective	Contact Heidolph Sales or Service
Shaking movement stops	Motor overheat protection has tripped	Wait about 20 minutes and reduce the load on the shaking platform
	Electrical fault (no engine noise)	Contact Heidolph Sales or Service
	Mechanical defect (engine noise audible)	Contact Heidolph Sales or Service

In case of recurring errors please contact the responsible sales department or our technical service (see section "Contact details" on page 46).

Technical specifications

Model-Specific data

HEI-MIX Titramax 100/101

Dimensions (W × H × D)	245 × 125 × 310 mm
Weight	5.5 kg
Permissible load	max. 2 kg
Usable area	220 × 220 mm
Motion	orbital vibrating
Orbit	Titramax 100: 1.5 mm Titramax 101: 3 mm
Speed range	150 – 1,350 rpm

HEI-MIX Titramax 1000

Dimensions (W × H × D)	320 × 128 × 375 mm
Weight	8 kg
Permissible load	max. 5 kg
Usable area	290 × 258 mm
Motion	orbital vibrating
Orbit	1.5 mm
Speed range	150 – 1,350 rpm

HEI-MIX Vibramax 100

Dimensions (W × H × D)	245 × 125 × 310 mm
Weight	5.5 kg
Permissible load	max. 2 kg
Usable area	220 × 220 mm
Motion	orbital vibrating
Orbit	3 mm
Speed range	150 – 1,350 rpm

HEI-MIX Rotamax 120


Dimensions (W × H × D)	245 × 125 × 310 mm
Weight	5.5 kg
Permissible load	max. 2 kg
Usable area	220 × 220 mm
Motion	orbital
Orbit	20 mm
Speed range	20 – 300 rpm

General device data**Hei-MIX Titramax 100/101/1000, Vibramax 100, Rotamax 120**


Drive	Shaded pole motor	
Speed control	analogue	
Protection class	IP30	
Degree of pollution	2	
Acoustic pressure	< 50 dB (A)	
Engine overheat protection	self-resetting	
Electrical data		
Rated voltage (variant-dependent)	Titramax 100/101	
	1/N/PE 230 V AC, 50/60 Hz	
	1/N/PE 115 V AC, 50/60 Hz	
	Titramax 1000	
	1/N/PE 230 V AC, 50 Hz	
Fuse	1/N/PE 115 V AC, 60 Hz	
	Vibramax 100 & Rotamax 120	
	1/N/PE 230 V AC, 50/60 Hz	
	1/N/PE 115 V AC, 50/60 Hz	
	115 V, 60 Hz	
Overvoltage category	M 1,25 A (230 V) or T 2,0 A (115 V)	
Protection class	II	
Permissible mains voltage fluctuations	I ⚡	
Power input	±10 %	
	Titramax 100/101	25 W
	Titramax 1000	25 W
	Vibramax 100	31 W
	Rotamax 120	33 W
Permissible ambient conditions		
Operating temperature	5 °C – 31 °C at up to 80 % rel. humidity	
	32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)	
Installation altitude	up to 2,000 m asl	

Performance range/maximum load


Hei-MIX Titramax 100/101, Hei-MIX Vibramax 100

HEI-MIX Titramax 1000

HEI-MIX Rotamax 120

Scope of delivery

Item	Quantity	Product no.
HEI-MIX Titramax 100*	1	544-11200-00
HEI-MIX Titramax 101*	1	544-11300-00
HEI-MIX Titramax 1000*	1	544-12200-00
Hei-MIX Vibramax 100*	1	544-21200-00
HEI-MIX Rotamax 120*	1	544-41200- 04
Power supply cord	1	country specific
Guarantee registration / Certificate of decontamination	1	01-006-002-78

* the product numbers shown are for 230 V devices for the European market. For product numbers of variants, please contact Heidolph instruments.

Accessories

HEI-MIX Titramax 100/101, Vibramax 100, Rotamax 120

	Quantity	Product no.
Device fuse for 230 V version	1	14-002-015-23
Device fuse for 115 V version	1	14-002-015-45

Accessories Hei-MIX Vibramax 100, Rotamax 120

Attachment with tension rollers	1	549-810000-00
Tension roller	1	11-008-007-08
Tablar 100 with universal perforation for use with clamps for Erlenmeyer flasks	1	549-59100-00

Further information on available accessories can be found on our website at www.heidolph.com.

Device service

When carrying out service work on the device (cleaning, maintenance, repair), observe the general instructions and safety information described in this section.

DANGER

Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

When cleaning, avoid the penetration of liquids.

Before replacing the fuses, switch off the device and disconnect the power supply cord.

Always replace the two device fuses in pairs with original manufacturer fuses. Further information on available accessories can be found on our website at www.heidolph.com.

After fuse replacement, check the device for a safe condition according to IEC 61010-1.

General cleaning instructions

Wipe all surfaces and the control panel with a damp cloth if necessary. Persistent contamination can be removed with mild soapy water.

CAUTION

Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

Never use any aggressive or abrasive cleaning agents or aids.

Replace the device fuse

The fuse holder for the two device fuses is located on the rear of the device below the IEC appliance inlet.

DANGER

Electric shock

Before replacing the fuses, switch off the device and disconnect the power supply cord.

Always replace the two device fuses in pairs with original manufacturer fuses, see section "Accessories" on page 43.

After fuse replacement, check the device for a safe condition according to IEC 61010-1.

Repairs - Return of equipment

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim. Regardless of the warranty claim, the owner is solely liable for damage caused by unauthorized repairs.

- In case of repair and before returning your device, contact our technical service at the following e-mail address:
 - service@heidolph.de.
- In your message, please provide us with the following information in addition to an error description:
 - Item number
 - Serial no.

The required data can be found on the rating plate of the device.

A service representative will contact you as soon as possible to agree on the next steps.

Include the completed certificate of decontamination with every device return, see "Certificate of decontamination" on page 47.

Maintenance

The device contains no user-serviceable components. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example) contact our technical service.

Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
- When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
- Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
- Properly remove and dispose of residues of substances that are hazardous to health, the environment and biohazardous!

Registration in accordance with the Electrical and Electronic Equipment Act (ElektroG)

Heidolph Scientific Products GmbH, headquartered in 91126 Schwabach/Germany, address Walpersdorfer Str. 12, is registered under the number DE 50705753 in the electronic waste equipment register of the foundation ear (www.ear-system.de).

Warranty statement

Heidolph Scientific Products GmbH provides a three-year warranty against material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage www.heidolph.com.

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

Contact details

Heidolph Scientific Products GmbH

Technical service

Walpersdorfer Str. 12

D-91126 Schwabach/Germany

Email: service@heidolph.de

Representations

To find your local Heidolph distributor please visit www.heidolph.com

Certificate of decontamination

Enclose the certificate of decontamination, duly completed, with your device return. Submissions without a certificate of decontamination cannot be processed!

CERTIFICATE OF DECONTAMINATION

IN CASE OF RETURNS

Please fill in the required fields.

Note: The sender must package the goods properly and appropriately for transport.

Heidolph Scientific Products GmbH

Walpersdorfer Straße 12
91126 Schwabach

Phone: +49 (0) 9122 9920-380

Fax: +49 (0) 9122 9920-19

E-Mail: service@heidolph.de

SENDER

Name _____ First name _____
 Company/institution _____ Department _____
 _____ Workgroup _____
 Address _____
 ZC/City _____
 Country _____ Phone _____
 Email _____

DEVICE DETAILS

Article number _____ Serial no. _____

Ticket number _____

Reason for sending in _____

Has the device been cleaned, decontaminated/disinfected? **Yes** **No** (Please mark as applicable)

If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous? **Yes** **No** (Please mark as applicable)

If yes, with which substances did the device come into contact?

LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product (from serial number on 200066855) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

Reax 2	541-21001-XX		
Reax top	541-10000-XX		
Reax control	541-11000-XX		
Reax 20/4	541-20004-XX		
Reax 20/8	541-20008-XX		
Reax 20/12	541-20012-XX		
Titramax 100	544-11100-XX, 544-11200-XX		
Duomax 1030	543-32105-XX, 543-32205-XX, 543-32210-XX		
Polymax 1040	543-42105-XX, 543-42110-XX, 543-42205-XX, 543-42210-XX		
Polymax 2040	542-40005-XX, 542-40010-XX		
Multi Reax	545-10000-XX		
Vibramax 100	544-21200-XX	Titramax 101	544-11300-XX
Rotamax 120	544-41200-XX	Vibramax 110	544-31200-XX
Titramax 1000	544-12200-XX	Inkubator 1000	549-90010-XX
Unimax 1010	543-12310-XX	Promax 1020	543-22332-XX
Unimax 2010	542-10020-XX	Promax 2020	542-20020-XX
Hei SHAKE Orbital Core	546-11XXX-XX		

Maschinenrichtlinie / Machinery Directive 2006/42/EG
EMV-Richtlinie / Electromagnetic Compatibility Directive 2014/30/EU
Delegierte (RoHS-) Richtlinie / Delegated (RoHS) Directive 2015/863/EU
Angewandte (harmonisierte) Normen / (Harmonized) Standards applied:
EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

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Jörg Ziel
Quality Manager

Declaration of Conformity

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,

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Description:

Shaker & Mixer with different working directions

The object of the declaration described above is in conformity with the relevant UK Statutory Instruments (and their amendments):

2008 No. 1597

The Supply of Machinery (Safety) Regulations

2008 2016 No. 1091

The Electromagnetic Compatibility Regulations 2016

2012 No. 3032

The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

and complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only)

ProductIP (UK) Ltd.

8. Northumberland Av.

London WC2N 5BY

Signed for and on behalf of Heidolph Scientific Products GmbH

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices and components by Heidolph Scientific Products GmbH fulfill the requirements of the EU directive 2011/65/EU (RoHS) and the extension 2015/863. All devices are compatible with the requirements of the RoHS Directive for the following materials:

Max. 0,01% of weight in homogene material for cadmium and max. 0,1% of the weight in homogeneous material for Lead, Mercury, Hexavalent chromium, Polybrominated biphenyls (PBB), Polybrominated diphenyl ethers (PBDE), Bis(2-ethylhexyl) phthalate (DEHP), Butyl benzyl phthalate (BBP), Dibutyl phthalate (DBP), Diisobutyl phthalate (DIBP).

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 04. September 2024

Handwritten signature of Andreas Hahn in blue ink.

Andreas Hahn
Managing Director

Handwritten signature of Jörg Ziel in blue ink.

Jörg Ziel
Quality Manager

China RoHS DECLARATION OF CONFORMITY

Heidolph Scientific Products GmbH has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Scientific Products GmbH products is 25 years.

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This table is created according to SJ/T 11364-2014

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
- O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
- X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Scientific Products GmbH may enter into further devices or can be used together with other appliances .

With these products and appliances in particular, Heidolph Scientific Products GmbH will not take responsibility for the EFUP of those products and appliances.

Schwabach, 04.September.2024

Andreas Hahn
Managing Director

Jörg Ziel
Quality Manager

© Heidolph Scientific Products GmbH

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Scientific Products GmbH PRODUCTS