Betriebsanleitung
Operating instructions
Notice d'instructions

Hei-SHAKE Unimax 1010

 **heidolph**
research made easy

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Seite 4 – 25

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Datum

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Firmenstempel

Translation of the original instructions
Page 30 – 51

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Contents

Introduction

| | |
|---|----|
| About this document..... | 30 |
| Typographic conventions..... | 30 |
| Copyright protection..... | 30 |
| Safety instructions in the official languages of the European Union | 30 |

General notes

| | |
|---|----|
| Directives applied, product certification..... | 31 |
| Copyright notice | 31 |
| California Residents..... | 31 |
| Residual risk | 31 |
| Intended use..... | 31 |
| Compliant use..... | 32 |
| Reasonably foreseeable misuse | 32 |
| Transportation..... | 32 |
| Storage..... | 32 |
| Acclimatization | 32 |
| Permissible ambient conditions | 32 |
| Energy saving potential and cost efficiency | 32 |

Safety

| | |
|---|----|
| General safety instructions..... | 33 |
| Electrical safety..... | 33 |
| Interfaces..... | 33 |
| Data security | 33 |
| Operational safety | 34 |
| Occupational safety | 34 |
| Electrostatic charging of liquid or powdered media..... | 34 |
| Personal protective equipment (PPE)..... | 34 |
| Environmental protection..... | 35 |
| Biohazard..... | 35 |
| Other regulations..... | 35 |

Device description

| | |
|---|----|
| Mechanical design, control panel, terminals | 36 |
|---|----|

Commissioning

| | |
|------------------------------|----|
| Positioning the device | 37 |
| Power supply | 37 |
| Attachments..... | 38 |

Operation

| | |
|---|----|
| Switching the device on/off | 39 |
| Automatic restart | 39 |
| Activating the [Automatic restart] function | 39 |
| Operating behavior with automatic restart | 39 |
| Deactivating the [Automatic restart] function..... | 40 |
| Operating modes..... | 40 |
| Continuous operation | 40 |
| Timer mode..... | 41 |
| Adjusting the shaking frequency in timer mode | 41 |
| Remote mode..... | 42 |
| Activating/deactivating remote access | 42 |
| Device activation via PC | 42 |
| Device activation via server-based software..... | 42 |

Troubleshooting

| | |
|----------------------|----|
| Troubleshooting..... | 43 |
| Error codes | 43 |

Appendix

| | |
|--|----|
| Technical specifications | 44 |
| Performance range/maximum load | 44 |
| RS232 (ASCII) interface commands | 45 |
| Scope of delivery..... | 47 |
| Device service..... | 48 |
| General cleaning instructions..... | 48 |
| Replace the device fuse..... | 48 |
| Repairs - Return of equipment..... | 49 |
| Maintenance..... | 49 |
| Disposal..... | 49 |
| Warranty statement..... | 50 |
| Contact details | 50 |
| Certificate of decontamination..... | 51 |

About this document

This operating manual describes the features and operation of platform shakers, Hei-MIX Unimax 1010 type from construction year 2023. The operating instructions are an integral part of the described device!

Typographic conventions

This document uses the following symbols, signal words, and highlights:

| | Description |
|---|---|
|  | <p>Warning symbols in combination with a signal word indicate dangers:</p> <p>DANGER Indicates an immediate dangerous situation. Failure to respect the indications will result in death or serious injury.</p> <p>WARNING Indicates a potential danger. Failure to respect the indications will result in serious injuries.</p> <p>CAUTION Indicates a potential hazard. If not avoided, damage to property and minor to moderate injuries can occur.</p> |
|  | <p>Mandatory signs are used to indicate important and useful information on handling a product. This information is used to ensure operational safety and to maintain the value of the product.</p> |
| [...] | <p>Square brackets identify labels of control elements on the device as well as labels and entries in software masks and user interfaces.</p> |
| → | <p>The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.</p> |

Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, as well as utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Scientific Products GmbH. Any violation is subject to compensation for damage.

Safety instructions in the official languages of the European Union

A summary of all safety instructions in the official languages of the European Union can be found in our Safety Guide for the product group Hei-SHAKE (shakers and mixers, Ref. 01-005-006-99). This document is available for download on our homepage in the most up-to-date version.

Directives applied, product certification

CE marking

The device meets all requirements of the following directives:

- European Machinery Directive 2006/42/EC
- EMC Directive 2014/30/EU

Copyright notice

The software implemented in this product is protected by copyright laws. The rights holder is Heidolph Scientific Products GmbH, Walpersdorfer Straße 12, 91126 Schwabach/Germany. Any open source software components in this software are excluded from our copyright. Further information is available in the service area on our website www.heidolph.com.

Violations of copyright (e.g. unauthorized use or modification of the software) may result in civil claims (e.g. omission, compensation for damages) and/or entail criminal penalties within the meaning of the legislation of the destination country.

California Residents

Important information for California residents regarding Prop 65. Please visit www.P65Warnings.ca.gov for more information.

Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

Intended use

The described platform shaker in this document has been specially developed for the following tasks:

- Shaking
- Mixing
- Emulsifying
- Suspending
- Separating
- Dissolving
- Staining

The areas of application of the described platform shaker include chemical, biological and environmental-analytical laboratory and research applications, applications of basic research and comparable facilities.

Due to its design, the device in its delivery condition may only be used in analytical processes or in laboratory-like conditions in the food, cosmetics, and pharmaceutical industries as well as other comparable industries that manufacture products intended for consumption by humans or animals, or for use on humans or animals.

Any other use of this device is not considered as intended!

Compliant use

The compliance of each individual application must always be evaluated by the user. Any additional measures necessary to ensure compliance are always the responsibility of the user.

Reasonably foreseeable misuse

Additional measures may be necessary for use under conditions or for purposes deviating from the intended use, and/or specific guidelines and safety regulations must be observed. Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is the sole responsibility of the operator.

All risks resulting from improper use are borne solely by the operator.

The device may only be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the sole responsibility of the operator!

Transportation

During transportation, avoid severe shocks and mechanical stresses that can cause damage to the device.

Keep the original packaging in a dry and protected place for later use.

Storage

Always store the product in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

Acclimatization

After each transportation and after storage under critical climatic conditions (e.g., high temperature difference between inside and outside), allow the device to acclimatize at room temperature for at least two hours to prevent possible damage due to condensation before commissioning it in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections (power supply, tubing) only after the device has been acclimatized!

Permissible ambient conditions

The device is **NOT** suitable for outdoor use! The device is **NOT** suitable for use in potentially explosive areas!

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

Energy saving potential and cost efficiency

Even in standby mode, small amounts of energy are used to supply individual device functions. However, especially over longer periods of nonuse, this resource consumption adds up and causes avoidable costs.

Therefore, always switch off the device when not in use, especially for long periods of time (devices with main switch: Main switch in position [O]; devices without main switch: Unplug) to effectively avoid unnecessary energy consumption and the resulting costs.

General safety instructions

Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.

Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself and, where applicable, on connected devices or the supply connections.

If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.

Only use the device in accordance with the regulations for intended use.

Electrical safety

Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.

Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).

Always use the power supply cord provided with the device.

Prior to use, check that the device and the power supply cord are free of visible damage.

Make sure that the power plug can be reached directly at any time, in order to disconnect the device from the power supply at any time without delay.

Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Scientific Products GmbH.

Always disconnect the device from the power supply system before carrying out any maintenance, cleaning or repair work

Interfaces

Safely separate extra low voltage inputs and outputs rated 25 V AC or 60 V DC or above according to EN 61140, or use double or reinforced insulation according to EN 60730-1 or DIN 60950-1.

Use only shielded connection cables. Connect the shield to the connector housing. Always cover unused interfaces (front and back) to prevent liquids from entering the device.

Data security

The user is responsible for ensuring data security when transferring data between the described device and other devices.

Only use secure networks for the data transfer and avoid use of critical infrastructure.

Only use high-quality shielded data cables for the data transfer.

For data transfer via USB, an industrial standard USB hub should be preferably used to ensure the most stable connection possible.

Operational safety

Operate the device under a closed ventilated fume hood when working with potentially hazardous substances (see EN 14175 and DIN 12924).

Do not make any unauthorized changes or modifications to the device!

Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!

Rectify malfunctions or faults on the device immediately.

Switch off and disconnect the device from the power supply, preventing reconnection, if it is not possible to eliminate the malfunction or rectify the fault immediately.

Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

Occupational safety

Always use the prescribed personal protective equipment (PPE) such as protective clothing, safety goggles, protective gloves, safety shoes, etc.

Do not operate any other devices in the immediate vicinity of the device ...

- which can generate electromagnetic fields in the frequency range between 9×10^3 Hz to 3×10^{11} Hz,
- which generate emission or radiation sources in the frequency range 3×10^{11} Hz to 3×10^{15} Hz (in the optical spectral range wavelengths from 1,000 μm to 0.1 μm),
- which generate ultrasonic or ionizing waves.

Do not process any substances that could release energy in an uncontrolled manner (e.g. self-ignition).

Do not process substances in which the energy input through mixing poses a danger.

Wipe off any liquid that may have spilled on the device immediately.

Electrostatic charging of liquid or powdered media

Liquid or powdered media as well as the laboratory vessels used can become dangerously electrostatically charged due to the shaking movement or when emptying or cleaning the vessel. The intensity of the charge generated depends on the properties of the medium, the flow rate, the process parameters as well as the size, geometry and material of the laboratory vessels (source: TRGS 727, Section 4).

The evaluation of a potential explosion hazard due to electrostatic charging and the definition of suitable safety measures are always the responsibility of the operator and are always based on the substances to be processed and the individual application.

Personal protective equipment (PPE)

The operating company must determine and provide the necessary PPE, depending on the respective application and the media and chemicals used.

The required instruction of the personnel is solely within the operating company's responsibility.

Environmental protection

When processing environmentally hazardous substances, take appropriate measures to avoid hazards to the environment.

The evaluation of corresponding measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking of the device with the biohazard warning symbol.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Scientific Products GmbH.

The operating company is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Mechanical design, control panel, terminals

Designation

Function

- | | |
|--|---|
| <p>1 Shaking platform</p> <p>2 Display</p> <p>3 [Start/Stop] button</p> <p>4 Mode button [Timer]</p> <p>5 Mode button [Speed]</p> <p>6 Mode button [Power Return]</p> <p>7 Push-and-rotary control</p> <p>8 RS232 interface (Sub-D, 9-pin) for activating the device in remote mode.</p> <p>9 Appliance inlet with fuse holder</p> | <p>Surface for placing laboratory vessels; different variants can be selected, see scope of delivery and accessories.</p> <p>Display of the operating parameters and messages.</p> <p>Pushbutton to activate/deactivate the shaking function. The button LED lights up when the shaking function is activated.</p> <p>Pushbutton to activate the [Program timer] function: After pressing the button, the button LED lights up. In this operating state, the device timer can be programmed with the push-and-rotary control in a range from 1 min – 99 h 59 min.</p> <p>Pushbutton to display/adjust the speed (shaking frequency) in Timer mode.</p> <p>Pushbutton to activate/deactivate the [Automatic restart] function.</p> <ul style="list-style-type: none"> ▪ Rotary function: To adjust the operating parameters [Motor speed] or [Timer], see mode buttons [Speed] and [Timer]. ▪ Push function: In remote mode, the connection to the external control unit is interrupted by pushing the control, to use the device again via the integrated control elements. |
|--|---|

Positioning the device

CAUTION

Improper positioning of the device may result in direct and indirect damage to property due to falling and/or spilling of fluids!

The surface must be even, non-slip, proper and sufficiently stable.

The device must be freely accessible at all times.

Make sure that a sufficient safety distance to moving components can be maintained on all sides.

Observe the generally valid rules for the safe positioning of work equipment.

The professional mounting and positioning of the device including the supplied accessories is within the sole responsibility of the operator!

Heidolph Scientific Products GmbH shall under no circumstances be liable for any personal injury or damage to property resulting from failure to observe the instructions for the correct installation and positioning of the device.

Power supply

DANGER

Observe the instructions in section "General safety instructions" on page 33, in particular Electrical safety.

Power supply cord with locking system

Only use the supplied power supply cord with locking system!

Connect the cable coupling of the cord [2] to the IEC appliance inlet [1] on the back of the device. Make sure that the locking system engages [3] correctly.

Connect the power supply cord to a properly secured mains socket-outlet.

Disconnecting the power supply cord

Disconnect the power supply cord from the mains socket-outlet.

Press in the locking lever on the cable coupling and disconnect the power supply cord from the IEC appliance plug on the back of the device.

Attachments

Use the various optional attachments to fix your laboratory vessels. Please observe the relevant information and instructions in the appendix to this document.

At low shaking frequencies, vessels such as Petri dishes can also be placed directly on the non-slip contoured rubber mat.

Always place individual vessels in the center of the shaking platform.


Always distribute several vessels evenly on the shaking platform.

Switching the device on/off

Use the main switch to switch the device on and off.

During the start-up, the Heidolph logo appears for about two seconds on the device display, at the same time, all button LEDs light up.

As soon as the device is ready for use, the last set shaking frequency is shown [rpm] and the button LEDs turn off.

Automatic restart

When the [Automatic restart] function is activated, the device switches on automatically upon return of the power supply after a power interruption and the speed is augmented to the last set value.

WARNING

Before activating the function, take appropriate measures to make sure that, in case of a possible unmonitored restart, no hazard for persons or risk of damage to property can occur.

In the as-delivered condition, the [Automatic restart] function is deactivated and has to be activated or deactivated manually.

The function can be activated and deactivated at any time, i.e. without interrupting a running process.

The evaluation of potential risks by a possible unmonitored automatic restart as well as of appropriate security measures is the sole responsibility of the user!

Activating the [Automatic restart] function

Push the mode button [Power return] to activate the function.

The information [Power Return On] as well as the warning symbol for the automatic restart is shown on the display for about two seconds:

Operating behavior with automatic restart

When the power is restored after a power interruption, a warning signal is emitted for about one second, at the same time, the warning symbol for the automatic restart flashes for about three seconds.

After this warning time has expired, the motor speed is raised again to the previously set value!

The timer function has to be newly programmed and started manually after each power interruption.

Deactivating the [Automatic restart] function

Push the mode button [Power return] to deactivate the function.

The information [Power Return off] as well as the warning symbol for the automatic restart is shown on the display for about two seconds:

Operating modes

The device can be used in continuous operation or with a timer (timer mode). Regardless of the selected operating mode, observe the following specific safety instructions:

WARNING

- Whenever possible, use closed vessels for processing corrosive, toxic or biohazardous substances and seal them safely.
- Increase the speed only gradually to the desired shaking intensity, especially with open vessels, and observe the fluid movements.
- Please note the maximum lifting capacity of the platform of your device, see section "Performance range/maximum load" on page 44.
- Lower the shaking intensity or reduce the total load on the platform when the device starts to vibrate during operation.
- Always ensure that the vessels are evenly distributed on the platform.
- Pay particular attention to an even distribution of the weight load when vessels of different sizes and/or differently filled vessels are placed on the platform at the same time!
- Before switching on the device, make sure that all vessels are properly fixed on the platform.
- Use suitable attachments for the used vessels. Further information on available accessories can be found on our website at www.heidolph.com!
- Always use the necessary and appropriate personal protective equipment!

Continuous operation

In continuous operation, the shaking movement can be started and stopped as required via the button [Start/Stop] after switching on the device.

- Place the vessel(s) on the platform. Observe the information given in section "Attachments" on page 38.
- Switch on the device (see "Switching the device on/off" on page 39).
- Use the push-and-rotary control to set the required shaking frequency (motor speed):
 - The set value is shown on the device display.
 - The shaking frequency is immediately raised or reduced to the set value.

The shaking frequency can be adjusted at any time during operation. When setting and maintaining the shaking frequency, the numerical value is shown on the display in orange. As soon as the preselected shaking frequency is reached, the numerical value is shown on the display in white.

- Press the [Start/Stop] button to activate the shaking operation.
- Press the [Start/Stop] button again to deactivate the shaking operation

Timer mode

In timer mode, a specific period of time can be set for a continuous shaking movement.

- Place the vessel(s) on the platform. Observe the position information given in section "Attachments" on page 38.
- Switch on the device (see "Switching the device on/off" on page 39).
- Press the [Timer] button briefly. As soon as the button LED lights up, a timer value in the range from 1 min – 99 h 59 min can be defined by means of the push-and-rotary control.
- Press the push-and-rotary control briefly to confirm and adopt the input value.
 - The set value and a hold sign (**II**) are shown on the display.

When setting the timer, the numerical value is shown on the display in orange.
After confirming the input value, the numerical value is shown on the display in white.

- Press the [Start/Stop] button to activate the shaking mode.
 - The timer starts automatically, the hold sign(**II**) disappears.
 - The LEDs of the [Start/Stop] and [Timer] buttons light up.
 - The remaining time until the expiration of the programmed time period is shown as follows:
 - up to a remaining time of one minute in hh:mm;
 - from a remaining time of one minute in mm:ss.
- During the running timer mode, the actual shaking frequency can be shown and, if necessary, adjusted at any time via the [Speed] button.

When the shaking mode stops, the timer is also stopped but not reset:

- A hold sign appears on the display (**II**).
- The timer continues to run when the shaking movement is restarted, starting from the last count value.

After the expiration of the timer, the device generates an acoustic signal (three short notification sounds) and the button LEDs [Start/Stop] and [Timer] turn off.

- Press and hold the [Timer] button for at least three seconds to deactivate the function:
 - The device can continue to be used in continuous operation.
 - The LED of the [Timer] button disappears.

Adjusting the shaking frequency in timer mode

In timer mode, the shaking frequency (motor speed) can be adjusted at any time via the [Speed] button:

- Press the [Speed] button during the running timer mode.
- The actual shaking frequency is shown on the display.
- Use the push-and-rotary control to set the required shaking frequency (motor speed):
 - The set value is shown on the device display.
 - The shaking frequency is immediately raised or reduced to the set value.

The shaking frequency can be adjusted at any time during operation. When setting and maintaining the shaking frequency, the numerical value is shown on the display in orange. As soon as the preselected shaking frequency is reached, the numerical value is shown on the display in white.

Remote mode

WARNING

Secure the device in remote mode with a clearly visible warning sign and, if necessary, take further ambient-specific protective measures that protect against damage to property and injuries in the event of unexpected/unmonitored starting of the device.

Activating/deactivating remote access

To activate the remote access of the device, a message of type "SetTargetSpeed", "SetMotorState", "SetTimerValue" oder "SetTimerState" has to be sent. The content of the first message is ignored. Afterwards, the remote symbol appears on the display and remote commands are executed normally.

The remote access can be stopped by pressing the push-and-rotary control.

Interface commands, see section "RS232 (ASCII) interface commands" on page 45.

Device activation via PC

The device can be addressed via the integrated RS232 interface via an external control. The data interface enables activation of the device and recording of the process data in remote mode. This requires suitable software such as Hei-PROCESSING Solutions of Heidolph Scientific Products GmbH, which supports the set of commands stored in the device.

Device activation via server-based software

The device can be addressed via the integrated RS232 interface via a server-based software. If required, contact an authorized dealer or our technical service, see section "Contact details" on page 50.

In remote mode, all local control elements with the exception of the [Start/Stop] button and the push function of the push-and-rotary control are deactivated:

- By pressing the [Start/Stop] button, the shaking mode can be stopped in remote mode.
- By pressing the push-and-rotary control, the external control is deactivated.

In both cases, external commands are ignored.

Troubleshooting

| Failure | Cause | Possible remedy |
|---|--|---|
| Display remains dark after switching on | No mains voltage | Check if the on/off switch is in position ON |
| | Power supply cord defective | Check the power supply cord for damages/correct connection |
| | On/off switch in position OFF | Check building fuses |
| | Display defective | Check device fuses |
| LED of the [Start/Stop] button does not light up when the function is activated | Button LED defective | Contact Heidolph Sales or Service |
| LED of the [Start/Stop] button does not light up when the function is activated | Button LED defective | Contact Heidolph Sales or Service |
| Shaking movement stops | Motor overheat protection has tripped | Wait about 20 minutes and reduce the load on the shaking platform |
| | Electrical fault (no engine noise) | Contact Heidolph Sales or Service |
| | Mechanical defect (engine noise audible) | Contact Heidolph Sales or Service |

In case of recurring errors please contact the responsible sales department or our technical service (see section "Contact details" on page 50).

Error codes


| Code | Description, remedy |
|------|---|
| E11 | Motor driver error Switch off and on the device. |
| E12 | Speed deviation between actual speed and nominal speed. ≥ 5 rpm or 10 % for at least 30 s. Switch off and on the device. |
| E21 | Timeout remote control, no valid message for at least 10 s. Set command "Reset Communication". |
| E41 | CRC error: calculated and stored CRC do not match. Switch off and on the device. |

Technical specifications

Hei-SHAKE Unimax 1010

| | |
|-----------------------------------|--|
| Dimensions (W × H × D) | 320 × 128 × 375 mm |
| Weight | 9.25 kg |
| Interface(s) | 1 × RS232 (Sub-D, 9 pole) |
| Permissible load | max. 5 kg |
| Usable area (W × D) | 290 × 258 mm |
| Movement / orbit | orbital / 10 mm |
| Speed range | 30 – 500 rpm |
| Drive | EC motor |
| Protection class | IP30 (acc. to EN 60529) |
| Protection class | I Ⓢ |
| Acoustic pressure | < 50 dB (A), in accordance with IEC 61010 |
| Rated voltage | 100 – 240 V, 50/60 Hz |
| Fuse | M 1.25 A (250 V) |
| Power input | 50 W |
| Operating and storage temperature | 5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly) |
| Installation altitude | up to 2,000 m asl |

Performance range/maximum load

RS232 (ASCII) interface commands

GET VERSION

With [Get Version], the version information can be read out.

→ Request: **STATUS[] r\n**

→ Response:

- Device Version: **STATUS DV= ID,VERSION,NAME,INFO\r\n**
- Firmware Version: **STATUS FW= ID,VERSION,NAME,INFO\r\n**
- Hardware Version: **STATUS HW= ID,VERSION,NAME,INFO\r\n**
- Variant version: **STATUS VV= ID,VERSION,NAME,INFO\r\n**

GET DEVICE STATE

With [Get Device State], the device state can be checked.

→ Request: **IN_MODE_ALL \r\n**

→ Response:

- **IN PV ALL**
 - TS [TARGET SPEED]**; (read out motor target speed)
 - AS [ACTUAL SPEED]**; (effective motor speed)
 - RH [TIME REMAINING HOURS]**; (remaining hours)
 - RM [TIME REMAINING MINUTES]**; (remaining minutes)
 - SH [TIME SET HOURS]**; (set value hours)
 - SM [TIME SET MINUTES]**; (set value minutes)
 - DOX [DEVICE]**; (device state*)
 - EC [ERRORCODE]**; (error code)
 - EL1 [ERRORCODELOG0]**; (error history)
 - EL2 [ERRORCODELOG1]**; (error history)
 - EL3 [ERRORCODELOG2]**; (error history)\r\n
- (*device status: 0 = Init, 1 = PowerOn, 2 = AutoRestart, 3 = Idle, 4 = Runtime, 5 = RunPermanent, 6 = IdleRemote, 7 = RunRemoteSpeed, 8 = RunRemoteTime, 9 = Error)

SET TARGET SPEED

With [Set Target Speed], the target speed of the motor can be set. This command is only executed if the system is in remote mode.

→ Request: **OUT_SP_3 [TARGET SPEED]\r\n**

→ Response: **OUT_SP_3 [TARGET SPEED]\r\n**

SET MOTOR STATE

With [Set Motor State], the motor can be switched on and off. This command is only executed if the system is in remote mode.

→ Request: **2 2\r\n**

→ Response: **START 2\r\n**

(**Set Motor State Stop** is used to switch off the motor.)

SET TIMER VALUE

With [Set Timer Value], the timer can be set to a value. This command is only executed if the system is in remote mode. By sending the command with the time 0 or 00:00 AM during operation, the timer ends and the motor continues turning until it is manually deactivated.

→ Request: **OUT_SP_9 HH:MM\r\n**

→ Response: **OUT_SP_9 HH:MM\r\n**

SET TIMER STATE

With [Set Timer State], the motor can be started or stopped for the set time. This command is only executed if the system is in remote mode.

- Request:
 - **OUT_MODE_9** 1\X\n (start timer or motor)
 - **OUT_MODE_9** 0\X\n (Stop timer or motor)
- Response: **OUT_MODE_9**\X\n

SET AUTOMATIC RESTART

With [Set Automatic Restart], the automatic restart can be defined after switching on the device. This command is only executed if the system is in remote mode.

- Request: **OUT_MODE_2** [**AUTOMATICRESTART**] r\n
- Response: **OUT_MODE_2** [**AUTOMATICRESTART**]\X\n

(0 = automatic restart deactivated, 1 = automatic restart activated)

SET CONNECTION CONTROL

With [Set Connection Control], the option [Remote Timeout Monitoring] can be activated or deactivated. In case of activation, the device quits the remote mode if it does not receive any messages over a period of ten seconds.

- Request:
 - **CC_ON**\X\n (Activate Connection Control)
 - **CC_OFF**\X\n (Disable Connection Control)
- Response:
 - **CC_ON**\X\n (Connection Control enabled)
 - **CC_OFF**\X\n (Connection Control disabled)

RESET COMMUNICATION

With [Reset Communication], a Remote Timeout Error can be reset. Afterwards, a remote connection can be established again.

- Request: **RESET**\X\n
- Response: **RESET**\X\n

GET ERROR STATE

With [Get Device State], the actual error state can be checked.

- Request: **ERROR**\X\n
- Response:
 - **ERROR EL1** [**ERRORCODELOG0**] ;
 - ERROR EL2** [**ERRORCODELOG1**] ; **ERROR EL3** [**ERRORCODELOG2**] ;\R\n

SET MOTOR STATE STOP

With [Set Motor State Stop], the motor can be switched off. This command is only executed if the system is in remote mode.

- Request: **STOP_2**\X\n
- Response: **STOP_2**\X\n

The transmission of the messages is made via ASCII strings. The protocol is active after the start. By a defined ADIP "Connect" command, you can switch to the ADIP mode.

RS232 interface parameters

- 9600 Baud
- Parity: none
- Data: 8 Bit
- Stop: 1 Bit

Command termination

All ASCII strings have to be finished with `\r\n` (0x0D 0x0A):

- ...`\r\n`

For further information on interface handling, please contact our technical service, see "Kontakt Daten" on page 4.

Scope of delivery

| Item | Quantity | Product no. |
|-----------------------|----------|------------------|
| Hei-MIX Unimax 1010 | 1 | 543-12310-00-8 |
| Power supply cord | 1 | country specific |
| Warranty registration | 1 | 01-006-002-78 |

Further information on available accessories can be found on our website at www.heidolph.com!

Device service

When carrying out service work on the device (cleaning, maintenance, repair), observe the general instructions and safety information described in this section.

DANGER

Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

When cleaning, avoid the penetration of liquids.

Before replacing the fuses, switch off the device and disconnect the power supply cord.

Always replace the two device fuses in pairs with original manufacturer fuses. Further information on available accessories can be found on our website at www.heidolph.com.

After fuse replacement, check the device for a safe condition according to IEC 61010-1.

General cleaning instructions

Wipe all surfaces and the control panel with a damp cloth if necessary. Persistent contamination can be removed with mild soapy water.

CAUTION

Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

Never use any aggressive or abrasive cleaning agents or aids.

Replace the device fuse

The fuse holder for the two device fuses is located on the rear of the device below the IEC appliance inlet.

DANGER

Electric shock

Before replacing the fuses, switch off the device and disconnect the power supply cord.

Always replace the two device fuses in pairs with original manufacturer fuses. Further information on available accessories can be found on our website at www.heidolph.com.

After fuse replacement, check the device for a safe condition according to IEC 61010-1.

Repairs - Return of equipment

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim. Regardless of the warranty claim, the owner is solely liable for damage caused by unauthorized repairs.

- In case of repair and before returning your device, contact our technical service at the following e-mail address:
 - service@heidolph.de.
- In your message, please provide us with the following information in addition to an error description:
 - Item number
 - Serial no.

The required data can be found on the rating plate of the device.

A service representative will contact you as soon as possible to agree on the next steps.

Include the completed certificate of decontamination with every device return, see "Certificate of decontamination" on page 51.

Maintenance

The device contains no user-serviceable components. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example) contact our technical service.

Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
- When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
- Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
- Properly remove and dispose of residues of substances that are hazardous to health, the environment and biohazardous!

Registration in accordance with the Electrical and Electronic Equipment Act (ElektroG)

Heidolph Scientific Products GmbH, headquartered in 91126 Schwabach/Germany, address Walpersdorfer Str. 12, is registered under the number DE 50705753 in the electronic waste equipment register of the foundation ear (www.ear-system.de).

Warranty statement

Heidolph Scientific Products GmbH provides a three-year warranty against material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage www.heidolph.com.

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

Contact details

Heidolph Scientific Products GmbH

Technical service

Walpersdorfer Str. 12

D-91126 Schwabach/Germany

Email: service@heidolph.de

Representations

To find your local Heidolph distributor please visit www.heidolph.com

Certificate of decontamination

Enclose the certificate of decontamination, duly completed, with your device return. Submissions without a certificate of decontamination cannot be processed!

CERTIFICATE OF DECONTAMINATION

IN CASE OF RETURNS

Please fill in the required fields.

Note: The sender must package the goods properly and appropriately for transport.

Heidolph Scientific Products GmbH

Walpersdorfer Straße 12
91126 Schwabach

Phone: +49 (0) 9122 9920-380

Fax: +49 (0) 9122 9920-19

E-Mail: service@heidolph.de

SENDER

Name _____ First name _____
 Company/institution _____ Department _____
 _____ Workgroup _____
 Address _____
 ZC/City _____
 Country _____ Phone _____
 Email _____

DEVICE DETAILS

Article number _____ Serial no. _____

Ticket number _____

Reason for sending in _____

Has the device been cleaned, decontaminated/disinfected? Yes No (Please mark as applicable)

If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous? Yes No (Please mark as applicable)

If yes, with which substances did the device come into contact?

LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product (from serial number on 200066855) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

| | | | |
|------------------------|--|----------------|--------------|
| Reax 2 | 541-21001-XX | | |
| Reax top | 541-10000-XX | | |
| Reax control | 541-11000-XX | | |
| Reax 20/4 | 541-20004-XX | | |
| Reax 20/8 | 541-20008-XX | | |
| Reax 20/12 | 541-20012-XX | | |
| Titramax 100 | 544-11100-XX, 544-11200-XX | | |
| Duomax 1030 | 543-32105-XX, 543-32205-XX, 543-32210-XX | | |
| Polymax 1040 | 543-42105-XX, 543-42110-XX, 543-42205-XX, 543-42210-XX | | |
| Polymax 2040 | 542-40005-XX, 542-40010-XX | | |
| Multi Reax | 545-10000-XX | | |
| Vibramax 100 | 544-21200-XX | Titramax 101 | 544-11300-XX |
| Rotamax 120 | 544-41200-XX | Vibramax 110 | 544-31200-XX |
| Titramax 1000 | 544-12200-XX | Inkubator 1000 | 549-90010-XX |
| Unimax 1010 | 543-12310-XX | Promax 1020 | 543-22332-XX |
| Unimax 2010 | 542-10020-XX | Promax 2020 | 542-20020-XX |
| Hei SHAKE Orbital Core | 546-11XXX-XX | | |

Maschinenrichtlinie / Machinery Directive 2006/42/EG
EMV-Richtlinie / Electromagnetic Compatibility Directive 2014/30/EU
Delegierte (RoHS-) Richtlinie / Delegated (RoHS) Directive 2015/863/EU
Angewandte (harmonisierte) Normen / (Harmonized) Standards applied:
EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

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Jörg Ziel
Quality Manager

Declaration of Conformity

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,

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Hei SHAKE Orbital Core 546-11XXX-XX

Description:

Shaker & Mixer with different working directions

The object of the declaration described above is in conformity with the relevant UK Statutory Instruments (and their amendments):

2008 No. 1597

2008 2016 No. 1091

2012 No. 3032

The Supply of Machinery (Safety) Regulations

The Electromagnetic Compatibility Regulations 2016

The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

and complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only)

ProductIP (UK) Ltd.

8. Northumberland Av.

London WC2N 5BY

Signed for and on behalf of Heidolph Scientific Products GmbH

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices and components by Heidolph Scientific Products GmbH fulfill the requirements of the EU directive 2011/65/EU (RoHS) and the extension 2015/863. All devices are compatible with the requirements of the RoHS Directive for the following materials:

Max. 0,01% of weight in homogene material for cadmium and max. 0,1% of the weight in homogeneous material for Lead, Mercury, Hexavalent chromium, Polybrominated biphenyls (PBB), Polybrominated diphenyl ethers (PBDE), Bis(2-ethylhexyl) phthalate (DEHP), Butyl benzyl phthalate (BBP), Dibutyl phthalate (DBP), Diisobutyl phthalate (DIBP).

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 04. September 2024

Handwritten signature of Andreas Hahn in blue ink.

Andreas Hahn
Managing Director

Handwritten signature of Jörg Ziel in blue ink.

Jörg Ziel
Quality Manager

China RoHS DECLARATION OF CONFORMITY

Heidolph Scientific Products GmbH has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Scientific Products GmbH products is 25 years.

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This table is created according to SJ/T 11364-2014

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Scientific Products GmbH may enter into further devices or can be used together with other appliances .

With these products and appliances in particular, Heidolph Scientific Products GmbH will not take responsibility for the EFUP of those products and appliances.

Schwabach, 04.September.2024

Andreas Hahn
Managing Director

Jörg Ziel
Quality Manager

© Heidolph Scientific Products GmbH

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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