
Fig. 1


Fig. 2

IKA

designed for scientists

IKA-Werke GmbH & Co. KG

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Phone: +84 28 38202142

eMail: sales.lab-vietnam@ika.com

Discover and order the fascinating products of IKA online:

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
Mounting instructions - IKA cross sleeves

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
R 271 Kreuzmuffe
R 271 Cross sleeve
R 271 Noix de serrage

20000004927_20008045a_Cross sleeve_082020

Technical specifications may be changed without prior notice.

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Safety instructions

- Read the mounting instructions in full before starting up and follow the safety instructions.
- Keep the operating instructions in a place where they can be accessed
- Safe operation of **IKA** laboratory devices with stands and cross sleeves is only guaranteed for combinations listed in the corresponding device operating instructions.
- Check that the laboratory device and the cross sleeve are properly tightened before each use.
- The laboratory device may only be moved on the stand when the device is at a standstill and unplugged!
- Clamped joints may loosen, resulting in danger to users and trial set-up!
- Hold the cross sleeve while loosening the clamped joints.
- While opening the clamped joints no sensitive parts are allowed to be under the device.
- Do not substitute turning handles or tighten with a tool.
- Do not use faulty cross sleeves!
- Observe the technical data!

Assembly

- Ensure that the extension arm on the laboratory device is tight; check regularly and retighten if necessary.
- Attach cross sleeve (**K**, see **Fig. 1**) to stand (**I**, see **Fig. 1**); open side of cross sleeve must face up (**R 182 / R 270**)!
- Mount laboratory device extension arm (**L**, see **Fig. 1**).
- When laboratory device is in desired position, firmly finger-tighten both turning handles (**M**, see **Fig. 1**) one after another.

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