[Gebruiksaanleitung](#) | [Operating manual](#) | [Mode d'emploi](#) | [Instrucciones de manejo](#) | [Istruzione](#) | [Instruções de utilização](#) | [操作手册](#) | [Руководство по эксплуатации](#) | [사용 지침](#) | [Bruksanvisning](#) | [使用説明書](#) | [Használati utasítás](#) | [Návod k použití](#) | [Gebruiksaanwijzing](#) | [Instrukcja użytkowania](#) | [Kullanım Talimatları](#)

HandyStep® S

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Otto-Schott-Str. 25
97877 Wertheim (Germany)

T +49 9342 808 0
F +49 9342 808 98000
info@brand.de
www.brand.de

Do you need more operating manuals and translations?
Please refer to <http://www.brand.de/om> or use the following
QR Code, then enter the product name and choose 'operating
manual' for document type.

The original operating manual is written in German. Other
languages are translations of the original operating manual.
Technical changes, errors or misprints reserved.

HandyStep® S is made in Germany under BRAND patents EP2415525B1,
EP2428271B1 and US8,652,418B2

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Contents

1	Using the operating manual	37
1.1	Signal words and their meaning	37
1.2	Symbols used in the instructions	37
1.3	Presentation of descriptions of actions	37
2	Scope of delivery	39
3	Safety Rules	40
3.1	Safety regulations	40
3.2	Function	41
3.3	Limitations of use	41
3.4	Operating exclusions	41
3.5	Limitations of use	41
4	Functional and operating elements	42
5	Operation	44
5.1	Inserting a PD-Tip	44
5.2	Setting the volume	46
5.3	Filling a PD-Tip	48
5.4	Repetitive dispensing	49
5.5	Ejecting a PD-Tip	50
5.6	Adapter for 25 ml and 50 ml	50
6	Checking the Volume (Calibration)	51
7	Cleaning and disinfection	53
7.1	UV sterilization	53
8	Troubleshooting	54
9	Product markings	55
10	Technical data	56
10.1	Accuracy table	56

10.2 Dispenser tips from other manufacturers	57
11 Order Information.....	58
12 Repairs.....	60
12.1 Sending for repair	60
13 Calibration service	62
14 Warranty	63
15 Disposal	64

1 Using the operating manual

- Please carefully read the operating manual before using the device for the first time.
- Keep the operating manual in an easily accessible place. It is part of the instrument.
- Be sure to include the operating manual if you transfer possession of this device to a third party.
- You can find updated versions of the instructions for use at www.brand.de/om.

1.1 Signal words and their meaning

⚠WARNING! ...

The signal word **WARNING** indicates a dangerous situation that, if not avoided, could result in death or serious injury.

⚠CAUTION! ...

The signal word **CAUTION** indicates a hazardous situation that, if not avoided, could result in moderate or minor injury.

NOTICE ...

The signal word **NOTICE** is used to address actions that are not related to physical injury. Example: Possible property damage.

The note shows operating tricks or other useful information

1.2 Symbols used in the instructions

This symbol warns of a potential hazard.

This symbol warns of hazardous media and the hazards that can arise from liquid splashes.

The symbol indicates that the instruction manual and safety information must be observed.

The symbol suggests that a face shield be worn during work.

The symbol suggests that eye protection be worn.

The symbol suggests that you to wear protective clothing.

The symbol suggests that protective gloves be worn.

1.3 Presentation of descriptions of actions

Format	Meaning
1. Task	Indicates a task.
a., b., c.	Indicates the individual steps of a task.

Format	Meaning
>	Indicates a prerequisite for a task.
⇒	Indicates a result of a completed task.

2 Scope of delivery

HandyStep® S, DE-M marked, with serial number, quality certificate, and shelf mount.

3 Safety Rules

3.1 Safety regulations

Please read carefully!

The HandyStep® S can be used in combination with hazardous materials, work processes and equipment. However, the operating manual cannot cover all of the safety issues that may occur in doing so. It is the user's responsibility to ensure compliance with the safety and health regulations and to specify the corresponding restrictions before use.

1. Every user must read and observe this operating manual before using the instrument.
2. Follow the general hazard instructions and safety regulations, e.g. wear protective clothing, eye protection and protective gloves. When working with infectious or hazardous samples, the standard laboratory rules and precautions must be adhered to.
3. Follow the instructions given by the reagent manufacturer.
4. Use the instrument only for dispensing liquids within the defined limitations and restrictions of use. Comply with the operating exclusions; see Operating exclusions, p. 41 . If in doubt, contact the manufacturer or supplier.
5. Always perform work in a manner that does not endanger users or other people. Avoid splattering. Use only suitable vessels.
6. Avoid touching the tip opening when working with aggressive media.
7. Never use force.
8. Do not make any technical modifications. Disassembly of the instrument is not permitted.
9. Always check that the device is in proper working condition before use. If an instrument fault is signaled, stop dosing immediately and follow the instructions in the chapter Troubleshooting, p. 54 . Contact the manufacturer if needed.

3.2 Function

HandyStep® S is a repetitive pipette for fast and simple repetitive liquid dispensing. The combination of PD-Tip size and the stroke setting made on the stroke setting dial enables volumes between 2 µl and 5 ml to be dispensed with maximum precision and accuracy.

3.3 Limitations of use

See Limitations of use, p. 41 .

3.4 Operating exclusions

- When the instrument is operated properly, the dosed sample comes in contact only with the tip and not with the HandyStep® S itself.
- The user is responsible for checking the suitability of the instrument for the intended use. Before using third-party dispenser tips, the user must verify their suitability and carry out a function/volume test!
- Do not use the instrument for dispensing liquids that corrode polypropylene (including cylinder PD-Tips), polyethylene (piston PD-Tips), LCP (piston 0.1 ml PD-Tip) or PC/PBT and polycarbonate.

NOTICE Instrument and tips are not autoclavable

(sterile PD-Tips see Order Information, p. 58). PD-Tips are disposable products.


Standard disinfectants can be used to disinfect the instrument, as long as they do not corrode PC/PBT and PC.

3.5 Limitations of use

Operating range *)	15 °C ... 40 °C (59 °F ... 104 °F)
Vapor pressure	Up to 500 mbar
Viscosity	20 mPa s at 50 µl PD-Tip 260 mPa s at 5 µl PD-Tip 977 mPa s at 1.25 µl PD-Tip

*) Additional temperatures upon request

4 Functional and operating elements

Volume table

Tip size (ml)		
Settings	Tip size (ml)	Steps
0.1	0.5	1
1	2	10
2	10	20
3	15	30
4	20	40
5	25	50
6	30	60
7	35	70
8	40	80
9	45	90
10	50	100
11	55	110
12	60	120
13	65	130
14	70	140
15	75	150
16	80	160
17	85	170
18	90	180
19	95	190
20	100	200
21	105	210
22	110	220
23	115	230
24	120	240
25	125	250
26	130	260
27	135	270
28	140	280
29	145	290
30	150	300
31	155	310
32	160	320
33	165	330
34	170	340
35	175	350
36	180	360
37	185	370
38	190	380
39	195	390
40	200	400
41	205	410
42	210	420
43	215	430
44	220	440
45	225	450
46	230	460
47	235	470
48	240	480
49	245	490
50	250	500

The volume table on back of the instrument is two-sided. To change the table, push the clip up, turn the table around and push the clip back down again.

Shelf/rack mount


The shelf mount is comes with adhesive strips for installation. It can be used without the strip in the Transferpette® S 6-position table stand.

5 Operation

5.1 Inserting a PD-Tip


a.

NOTICE PD-Tips are disposable products!

- a. Push the locking/filling lever to the lower stop.

b.

- b. Swing out the locking/filling lever.


c.

- c. Insert PD-Tip straight and upright from the bottom until it clicks into place.

d.

- d. Swing back the locking/filling lever.

5.2 Setting the volume

The possible dispensing volumes are listed under the transparent clip on the back of the instrument. You can set the stroke setting to 9 positions. These positions determine how many steps can be dispensed per volume. The following table shows the relationship between stroke setting (Setting), nominal volume (Tip Size), dispensing steps (Steps), and step volume (Volume).

Setting	Tip size (ml)										Steps
	0.1	0.5	1	1.25	2.5	5	10	12.5	25	50	
1	2	10	20	25	50	100	200	250	500	1000	49
1.5	3	15	30	37.5	75	150	300	375	750	1500	32
2	4	20	40	50	100	200	400	500	1000	2000	24
2.5	5	25	50	62.5	125	250	500	625	1250	2500	19
3	6	30	60	75	150	300	600	750	1500	3000	15
3.5	7	35	70	87.5	175	350	700	875	1750	3500	13
4	8	40	80	100	200	400	800	1000	2000	4000	11
4.5	9	45	90	112.5	225	450	900	1125	2250	4500	10
5	10	50	100	125	250	500	1000	1250	2500	5000	9
Volume (µl)											

Examples: A dispensing volume of 200 µl is required per dispensing step.

PD-Tip size	Stroke adjustment	Dispensing step
2.5 ml	4	max. 11
5 ml	2	max. 24
10 ml	1	max. 49

Setting the volume

To adjust the volume, follow these steps:

- a. Search for the desired dispensing volume (volume) in the table.
- b. Use the stroke setting dial to set the stroke setting for the volume that corresponds to the PD-Tip (tip size).

Dispensing steps (Steps)

Depending on the PD-Tip used, the number of steps varies even when the dispensing volume is the same. The fewer the number of steps, the greater the accuracy of the dispensed volume.

The volumes listed in the table are the only options possible. Since the first step must be discarded, the number listed in the table is always 1 step less.

Special feature of the 25 ml and 50 ml PD-Tips

To achieve the tolerances for these PD-Tip sizes, proceed as follows:


- a. Fill the PD-Tip completely.
- b. Set the stroke setting to 3.
- c. Discard the 1st step.
- d. Then set the desired stroke setting.

Empty residual liquid

- a. Hold the PD-Tip tip at an angle against a vessel wall.
 - b. Push the locking/filling lever all the way down.
- ⇒ The PD-Tip empties.

5.3 Filling a PD-Tip


a.

NOTICE We recommend priming a new PD-Tip before use. To do this, fill the tip with a minimal amount of liquid and empty it again. Small air bubbles in the area of the piston after priming do not affect the results.

a. Push the locking/filling lever to the lower stop.

b.

b. Immerse the head of the PD-Tip vertically into the liquid 3 - 10 mm.

c.


c. Slowly pull the locking/filling lever up to the upper stop to prevent the formation of air bubbles.

NOTICE Small air bubble near the piston can be ignored since a locking mechanism prevents the residual volume from being dispensed inadvertently after the final dispensing step.

5.4 Repetitive dispensing


b.

NOTICE The first dispensing step must be discarded!

- a. Check the volume setting once again. (Compare with the settings in the volume table!)
- b. Wipe off liquid adhering to the outside of the PD-Tip tip with lint-free cellulose tissue.
- c. Place the PD-Tip tip against the vessel wall.

c.


d.


- d. Dispense the liquid by pressing down the dispensing lever completely, then let it slide back up again.
- e. Make sure to dispense smoothly and evenly.

5.5 Ejecting a PD-Tip


b.

c.

d.

⚠WARNING! Media contact possible

Media residues may be present on the PD-Tip tip. Only eject a PD-Tip in a manner that does not endanger the user or other people.

NOTICE If viscous media is repeatedly dispensed, the leak tightness of the PD-Tip can no longer be guaranteed.

- a. Hold the HandyStep® S over a vessel.
 - b. Empty the PD-Tip by pushing the locking/filling lever to its lower stop.
 - c. Swing out the locking/filling lever (this will cause the piston to decouple).
 - d. Hold the HandyStep® S over a waste container. Press the ejector button (releases the cylinder).
- ⇒ PD-Tip is ejected.

5.6 Adapter for 25 ml and 50 ml

Adapter

To use the 25 ml and 50 ml PD-Tips, the reusable adapter supplied with the tips is required. This is coupled to the PD-Tips with a bayonet lock and can be detached after use. Adapters are autoclavable at 121 °C (2 bar), in accordance with DIN EN 285.

6 Checking the Volume (Calibration)

Depending on use, we recommend that gravimetric testing of the instrument be carried out every 3-12 months. This time frame should be adjusted to correspond with individual requirements. Gravimetric volume testing according to DIN EN ISO 8655-7 is performed as follows:

1. Prepare instrument

Insert a tip. This test can be carried out with PD-Tips of any size. However, a 5 ml PD-Tip is normally used for this purpose.

2. Perform test

- a. Set stroke setting 5 on the HandyStep® S.
- b. Fill the PD-Tip. Immerse the PD-Tip vertically into the test liquid.
- c. Discard the first dispensing step (step). It serves to compensate for play.

NOTICE For 25 ml and 50 ml PD tips:

To achieve the tolerances, proceed as follows:

- Fill the PD-Tip.
 - Set stroke setting to 3 before discarding.
 - Discard the 1st step.
 - Set stroke setting to 5 and proceed as described.
- d. Place a weighing vessel filled with a small amount of deionized water on the scale and tare the scale.
 - e. Dispense the second step into the weighing vessel. In doing so, push the dispensing lever down at a uniform speed until the first stop and hold it. Then wipe the tip off over a distance of approx. 10 mm.
 - f. Enter the weight into the test record.

- g. Repeat points D to F a total of 10 times. (For the 10th repetition with the 5 ml PD-Tip, the PD-Tip must be refilled with stroke setting 5 – points A to C).
- h. Repeat the same testing procedure at stroke setting 3 (60% of nominal volume) and 1 (20% of nominal volume).
- ⇒ This results in a total of the 30 required weighed values.

Calculation (for nominal volume)

x_i = weighing results

n = number of weighings

V_0 = nominal volume

Z = Correction factor (e.g. 1.0029 µl/mg at 20°C, 1013 hPa)

Mean:

$$\bar{x} = \frac{\sum x_i}{n}$$

Mean volume:

$$\bar{V} = \frac{V_1 + V_2 + \dots + V_n}{n}$$

Accuracy*:

$$A\% = \frac{\bar{V} - V_0}{V_0} * 100$$

Standard deviation*:

$$s = Z * \sqrt{\frac{\sum (x_i - \bar{x})^2}{n - 1}}$$

Coefficient of variation*:

$$CV\% = \frac{100 s}{\bar{V}}$$

*) Accuracy and coefficient of variation are calculated according to the formulas of statistical quality control.

NOTICE Test Instructions (SOPs) are available for download from www.brand.de.

7 Cleaning and disinfection

7.1 UV sterilization

The device is resistant to normal exposure to a UV disinfection lamp. The device may change color due to the UV exposure. We recommend the following parameters for UV disinfection:

Light spectrum: UV-C

Wavelength: 220 nm ... 270 nm

Irradiation time per cm^2 and intended degree of inactivation: 2 s ... 300 s

8 Troubleshooting

Problem	Possible cause	Corrective action
PD-Tip piston cannot be locked after the PD-Tip has been inserted.	Filling/locking lever is not pushed down to bottom position and is not swung out.	First press the ejector button and remove the PD-Tip, then push the filling/locking lever all the way down and swing it out completely.
Undefined dispensing volume	The stroke setting dial is not properly engaged.	Ensure that the stroke setting dial is securely locked into place.
Air bubble under the PD-Tip piston increasing in size.	PD-Tip is leaking.	Replace PD-Tip.
A and/or CV exceed tolerances.	PD-Tip is leaking.	Replace PD-Tip.
	Dispenser tips are difficult to move.	Replace PD-Tip or use other more pliant dispenser tips.
	Dispensing lever has not been pressed completely and with even pressure.	Press the dispensing lever with even pressure to the stop position.
	The first dispensing step was not discarded.	Discard the first dispensing step.
PD-Tip cannot be decoupled after emptying completely.	The PD-Tip has been inserted incorrectly.	Press the ejector button and then push the filling/locking lever all the way down and swing it out.

9 Product markings

Symbol or number	Meaning
	General warning sign
	Note user manual
XXZXXXXX	Serial number
 25	<p>The device is marked in accordance with the German Weights and Measures Act and the Weights and Measures Ordinance.</p> <p>Character sequence DE-M (DE for Germany), framed by a rectangle, as well as the two last digits of the year the marking was added.</p>

10 Technical data

10.1 Accuracy table

PD-Tips // from BRAND, 20 °C 'Ex'

PD-Tips sizes	Volume range (µl)	A ≤ ± % *			*CV ≤ % *		
		The stroke setting corresponds to the percentage of the nominal volume**			Stroke setting corresponds to % of nominal volume**		
		1 ± 20 %	3 ± 60%	5 ± 100 %	1 ± 20 %	3 ± 60%	5 ± 100 %
0.1 ml	2...10	4	2.4	1.6	6	3	2
0.5 ml	10...50	2.5	1.5	1	2.5	1.5	1
1.0 ml	20 - 100	2.5	1.5	1	2	1.2	0.8
1.25 ml	25...125	2.5	1.4	0.9	2	1.1	0.7
2.5 ml	50...250	1.8	1.1	0.7	1.5	0.9	0.6
5.0 ml	100 - 500	1.8	1.1	0.7	3	1	0.6
10.0 ml	200 - 1000	1.8	1.1	0.7	2	0.66	0.4
12.5 ml	250...1250	1.8	1.1	0.8	2	0.66	0.4
25.0 ml	500 - 2500	1.5	0.9	0.6	1.5	0.5	0.3
50.0 ml	1000 - 5000	1.5	0.8	0.5	1.5	0.5	0.3

*A = Accuracy, CV = Coefficient of Variation

** The nominal volume refers to the maximum dispensable volume of the entire system, which consists of the multi-dispenser and the PD-Tip being used.

Error limits relative to the set dispensing volume as a function of PD-Tip size, at the same temperature (20 °C) of instrument, tip, ambient environment, and distilled water, as well as smooth and even handling. The test is carried out in accordance with DIN EN ISO 8655-7.

10.2 Dispenser tips from other manufacturers

We recommend using BRAND PD-Tips //; however, the HandyStep® S can also be used with compatible dispenser tips from other manufacturers. When using compatible dispenser tips, the user must verify their suitability and carry out a function/volume test!

11 Order Information

HandyStep® S

Description	Packaging unit	Item no.
HandyStep® S, 1 shelf mount	1	705110

Shelf/rack mount

Description	Packaging unit	Item no.
Shelf/rack mount	1	705130

PD-Tip II

Volume [ml]	Item no.	Packaging unit [pieces]	Item no. BIO-CERT® LIQUID HANDLING STERILE	Packaging unit [pieces]
0.1	705700	100	705730	100
0.5	705702	100	705732	100
1	705704	100	705734	100
1.25	705706	100	705736	100
2.5	705708	100	705738	100
5	705710	100	705740	100
10	705712	100	705742	100
12.5	705714	100	705744	100
25	705716	25+1 adapter	705746	25+1 adapter
50	705718	50+1 adapter	705748	25+1 adapter
Set PD tips II 0.5 ml ... 12.5 ml	705720	20 per set	—	—

BIO-CERT® LIQUID HANDLING QUALITY: Free from human DNA, RNase, DNase, pyrogens, PCR inhibitors and ATP.

Quality grades for tips at BRAND:

The BIO-CERT® quality grades are proprietary standards of BRAND GMBH + CO KG. They demonstrate compliance with our stringent internal product standards.

PD-Tip II Set

Description	Item no.
PD-Tip II Set, 20 PD-Tips II each in sizes 0.5, 1.0, 1.25, 2.5, 5.0, 10 and 12.5 ml	705720

Adapter for 25 ml and 50 ml PD-Tip II

Volume [ml]	Item no.	Packaging unit	Feature
25 ml and 50 ml	702398	10	
25 ml and 50 ml	702399	5	BIO-CERT® LIQUID HANDLING STERILE

BIO-CERT® LIQUID HANDLING QUALITY: Free from human DNA, RNase, DNase, pyrogens, PCR inhibitors and ATP

12 Repairs

12.1 Sending for repair

NOTICE Transporting hazardous materials without approval is prohibited by law.

Clean the instrument thoroughly and decontaminate!

- When returning products, please enclose a general description of the type of malfunction and the media used. If information regarding media used is missing, the instrument cannot be repaired.
- Shipment is at the risk and the cost of the sender.

Outside USA and Canada

Complete the “Declaration on Absence of Health Hazards” and send the instrument to the manufacturer or supplier. Ask your supplier or manufacturer for the form. The form can also be downloaded from www.brand.de.

Within USA and Canada

Please clarify the requirements for the return delivery with BrandTech Scientific, Inc **before** sending the instrument in for service.

Return only cleaned and decontaminated instruments to the address provided with the Return Authorization Number. Place the Return Authorization number so that it is clearly visible on the outside of the package.

Contact addresses

Germany:

BRAND GMBH + CO KG
Otto-Schott-Straße 25
97877 Wertheim (Germany)
T +49 9342 808 0

USA and Canada:

BrandTech® Scientific, Inc.
11 Bokum Road
Essex, CT 06426-1506 (USA)
T +1-860-767 2562

F +49 9342 808 98000
info@brand.de
www.brand.de

India:

BRAND Scientific Equipment Pvt. Ltd.
303, 3rd Floor, 'C' Wing, Delphi
Hiranandani Business Park,
Powai
Mumbai-400 076 (India)
T +91 22 42957790
F +91 22 42957791
info@brand.co.in
www.brand.co.in

F +1-860-767 2563
info@brandtech.com
www.brandtech.com

China:

BRAND (Shanghai) Trading Co., Ltd.
Rm 201-202, North Tower,
No. 199 Kaibin Rd, Xuhui District, Shanghai
Shanghai 200030 (P.R. China)
T +86 21 6422 2318
F +86 21 6422 2268
info@brand.com.cn
www.brand.cn.com

13 Calibration service

The ISO 9001 and GLP guidelines require regular inspection of your volume measuring devices. We recommend performing a volume check every 3 to 12 months. The cycle is dependent on the individual requirements of the instrument. Checks should be performed more frequently, in case of high frequency of use or the use of aggressive liquids.

The detailed testing instructions are available for download on www.brand.de or www.brandtech.com.

BRAND also offers the possibility to have your instruments calibrated by our factory calibration service or our accredited calibration laboratory. Simply send us the device to be calibrated, accompanied by details about which type of calibration you would like. The instruments will be returned to you after a few days. The instruments are supplied with a detailed calibration certificate or a calibration slip in accordance with DIN EN ISO/IEC 17025. More information can be obtained from your dealer or directly from BRAND. The order document is available for download on www.brand.de (see Service and Support).

For customers outside Germany

If you would like to use our calibration service, please contact one of our service partners in your region. If a factory calibration is requested, the service partner can forward the instrument to BRAND.

14 Warranty

We shall not be liable for the consequences of improper handling, use, servicing, operating or unauthorized repairs of the device or for the consequences of normal wear and tear, especially of wearing parts such as pistons, seals, valves and the breakage of glass. The same applies for failure to follow the instructions of the operating manual. We are not liable for damage resulting from disassembly beyond that described in the operating manual or if non-original spare parts or components have been installed.

USA and Canada:

Find more warranty information on www.brandtech.com.

15 Disposal

At the end of their service life, ensure the instruments are properly disposed of and, where possible, recycled. Disposal must be carried out in accordance with applicable national and local regulations.

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