S-Viewer

Mikroskop-Kamera-Software D

Seite 02 - 33

Microscope Camera Software GB

Page 34 – 66

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Operating Manual Camera Software

Table of Contents

1	Short Description.....	4
2	Features.....	4
3	Functions.....	5
4	Software Download.....	5
5	Run the software.....	5
6	Select the microscope/camera and object preview.....	6
6.1	Select the microscope/camera.....	6
6.2	Select the display size of the object preview	6
6.3	Select magnification.....	7
7	Operation of the software.....	8
7.1	Menu bar.....	8
7.2	File.....	8
7.3	Process.....	9
7.4	Methods of operation	10
7.5	View.....	12
7.6	Advanced	12
7.7	Count	13
7.8	Rotate.....	13
7.9	Stitching	13
7.10	High Dynamic Range (HDR)	16
7.11	Dynamic EDF	17
7.12	Mix.....	17
7.13	Automatic Stitching.....	18
7.14	Semi-automatic Stitching	18
7.15	Layer	19
7.16	Calibration.....	21
7.17	Show/Hide Scale Bar.....	22
7.18	Measure	23
7.19	Measurement settings.....	25
7.20	Export Data	25
7.21	Report	25
7.22	Annotation.....	27
7.23	Save settings	28
7.24	Language settings.....	29
7.25	Arrangement of open windows.....	29
7.26	Additional features	30
7.27	Interval capture.....	30
7.28	Video recording	31
8	Troubleshooting.....	32
8.1	Preview failed, blank screen.	32

1 Short Description

S-Viewer is a professional software specifically designed for digital microscope applications. It plays a central role in the interaction between the user and a digital microscope, meeting the requirements for visualizing operational structures on a PC.

2 Features

- Compatible with all USB 2.0 and USB 3.0 interfaces, as well as the KERN Optics microscope series OIV-3 / OIV-9 / OBT-2.
- Real-time video editing functions.
- Supports the measurement of dynamic videos, reduces the number of measurement steps, and shortens the measurement duration.
- With a customizable RGB adjustment that can be saved, it enables lifelike color reproduction. Additionally, multiple RGB options can be configured, enhancing user comfort and providing an optimal solution for color accuracy.
- Easy to use, it allows the embedding of measured data, file headers, dates, and scales into images, and supports quick storage.
- Ermöglicht das Hinzufügen detaillierter Dateiinformationen bei gleichzeitiger Bewahrung des Bildes, wodurch die Beschreibung von Bilddateien verbessert wird.

3 Functions

- Video preview, photography, continuous shooting
- Parameter settings such as brightness, contrast, saturation, chroma, gamma, sharpness
- Automatic white balance (AWB) and automatic exposure
- Real-time video processing such as horizontal and vertical mirroring, grayscale, embossing effect, sharpening, binary, negative image
- RGB three-color real-time adjustment, savable, controllable
- Full-screen mode
- Image counting, image stitching, HDR, EDF, rotation
- Performing a measurement calibration
- Dynamic real-time measurement
- Static image measurement

4 Software Download

The camera software can be downloaded for free from the KERN & SOHN website under Downloads (<https://www.kern-sohn.com/shop/de/downloads/>).

5 Run the software


After successfully downloading and installing the software, launch it via the icon

6 Select the microscope/camera and object preview

6.1 Select the microscope/camera


After connecting the microscope or camera to your computer via a USB cable, select the correct device from the "Choose Input Device" menu.

The connection will then be established, and you will see the object from the object table displayed on the screen.

6.2 Select the display size of the object preview

In the toolbar below the menu bar...


...or by holding down the Ctrl key and using the mouse scroll wheel, you can select the display size of the object preview.

...or switch to full-screen mode using this icon

6.3 Select magnification


Also, in the toolbar below the menu bar...

...switch between the saved magnification levels, see calibration.


7 Operation of the software

7.1 Menu bar

7.2 File

Under the "File" menu, various functions are available to manage images and control the application.

New Preview: Start a new object preview.

Open Image Folder: Opens the location where the images are stored.

Open...: Open saved images.

Close: Closes the current preview.

Save: Saves the image in the default folder.

Save As...: Allows you to choose and set storage attributes.

Print...: Available in static view mode to print the current image.

Print Preview: Displays a preview of the printed layout on A4 paper in static view mode.

Print Setup: Allows adjusting the print settings in static view mode.

Select Source...: Choose from the available cameras.


Acquire...: Take a picture in TWAIN mode.

Exit: Closes the application.

7.3 Process

Under this menu, you will find various real-time editing methods.


Pull-down menu

Tool bar

By activating or using the symbol with the "E," you can enable and disable the remaining functions.

7.4 Methods of operation

Original image


Horizontal flip


Vertical flip

Monochrom


Negative image


Sharpen


Original image


Embossing

Binary image

To obtain a dynamic binary representation of your sample, adjust the threshold for binarization (from 1 to 255). To do this, check the box and move the slider:


Start or stop the RGB real-time alternate view of the video by checking or unchecking the box here and adjusting the sliders:

The lower section will be displayed here when you click the "*Advanced*" button:

Apply: Select options from the dropdown list and apply them to the current video.

Add: Click "Add," enter the option names, and add the adjusted Red-Green-Blue values to the database. These will then appear in the list.

Delete: Select options from the dropdown list and delete them from the list or the database.


7.5 View

Through this menu, different areas on the display can be enabled or disabled. If the checkbox is checked, the function/information view will be displayed. If unchecked, it will be hidden. The pin allows you to pin the respective view window in the S-Viewer window , or be removed again .


7.6 Advanced

In this menu, you will find additional functions for editing your sample.


7.7 Count

The watershed algorithm, or watershed segmentation, is a technique used to segment images by separating regions within an image based on differences in brightness. The algorithm helps divide an image into different areas by identifying the "boundaries" between them.

Counting settings refer to the configurations that determine how objects or elements are counted in an image or analysis. These settings define the criteria that must be met for an object to be counted and how the counting process is carried out.

7.8 Rotate


This application can only be used with static images. The images can be rotated by 90 degrees, 180 degrees, 270 degrees, or an arbitrary angle clockwise.


7.9 Stitching

In this editing method, multiple individual images are combined into a larger, seamlessly connected image. The software automatically detects the overlapping areas and merges the images into a single, composite image.


Use  to manually create an image section. Now move the sample slightly and create another image. Repeat this process until you have photographed the entire sample area. All images are shown one after the other in the list (here IMAGE1, IMAGE2, etc.).


Stitching can be found via the menu item *Advanced...*


...and then under *Stitching*.


The images from the bar are now listed in the left-hand side of the window. Alternatively, you can also load saved images into the programme using the *Loading...* button.


Then press *Add All >>* and the listed images will appear on the right-hand side of the window.

Then press *OK* and the merging of the images will start. This process may take a few seconds.

At the end, the merged image (here *Stitching_0000.jpg*) is shown in the preview window and you can save the file to a location of your choice.


7.10 High Dynamic Range (HDR)

This is a technique that allows images with an extended brightness range to be created, so that both very bright and very dark details in an image are more visible.

Static HDR (Activate during dynamic image display) refers to a feature that allows the benefits of HDR (High Dynamic Range) to be utilized in a static image representation while displaying a dynamic (moving or live) image.

Disable the camera's automatic exposure, dim the light source to the appropriate brightness, click "*Start HDR Synthesis*", and then slowly adjust the brightness of the light source.

Dynamic HDR (Activated when displaying non-dynamic images):

Dynamic HDR refers to the application of HDR techniques while images are processed in real-time or in a dynamic display. This feature allows for optimizing brightness and details in images during live display, even with changing lighting


conditions or movement in the image. It is activated when the image is not static, meaning during moving or interactive displays.

Click this function to load the opened image by default into the left list. Click "*Loading*" to import an HDR image from an external source. Select the image that requires HDR, click "*Add*" to import it into the right list, and click "*OK*" to synthesize the selected image with HDR.

7.11 Dynamic EDF

Dynamic EDF (Activate during dynamic image display):

Refers to a feature that enhances image display by applying Dynamic Extended Depth of Field (EDF) during the display of dynamic (live or moving) images.


It can be opened in dynamic view and is available when automatic exposure and focus are turned off. The small image on the left is a real-time micro image, and the right image is a depth-of-field composite image. Adjust the focal length so that the image on the left is

sharp in a specific area and wait for the image on the right to be synthesized. After synthesis, adjust the focal length again to synthesize the remaining areas until all areas are clearly displayed. Once finished, click "Close" in the menu to obtain the EDF composite image.


7.12 Mix

You can blend two images according to different transparency ratios.

7.13 Automatic Stitching

Automatic stitching can be found under the *Advanced* menu item. Then select *Start* and you will see this overview with two images as shown here.


The left-hand image is the live image and the right-hand image is the composite result. Now move the object slowly in one direction and you will see on the right-hand image how it slowly assembles. If you move too quickly, the colour of the frame will change from green to red

and you will have to turn the live image back a little until the camera has found the correct position. The frame will turn green again. Repeat this process until you have created a complete image. Then select *Stop* in the menu under *Automatic stitches* and you can save the composite image.


7.14 Semi-automatic Stitching

You can also start semi-automatic stitching under the *Advanced* menu item. Select *Start* again and you will see this overview.


The rest of the procedure is very similar to the automatic method, with the difference that you trigger the images yourself using the space bar on your PC keyboard. When you have finished, select *Stop* again and you can save the entire image.

7.15 Layer


After starting the program and selecting the camera, the first image is set as the background layer by default. Additional layers can be added, deleted, and edited through the menu.

To add a new layer, click on "New..." and a new window will open to assign a name to this layer:

To delete a layer, select "Delete..." from the menu and click on the layer you want to delete:


To delete all layers, select "Delete ALL" from the menu.

Note: The background layer and the current layer cannot be deleted.


"Current...": Select the desired layer and set it as the current layer.

You can control the visibility of layers through the "*Show/Hide...*" menu. Check or uncheck the boxes to show or hide the layers.

Note: The background layer and the current layer cannot be hidden.



You can edit the current layer by selecting "*Rename...*"

7.16 Calibration

As preparation for later measurements, a calibration must be performed. This function can be found under the "*Measure*" menu. Clicking on it will open this window, and a green measurement line will appear on the display (right image).


Possible ruler units: μm , mm, cm, dm, inch, mil

Place a physical ruler (micrometer, e.g. ODC-A2404) on the objective stage and select the correct objective magnification (top left image, here 4x). Move the green measurement line (right image) and enter this value in the left image under "*Ruler Unit*" (here mm) and "*Actual Length*" (here 1 mm). Then confirm by clicking "OK." The values are now saved for this magnification. Repeat this process for all required objective magnifications.

If you want to replace a saved value, press the "*Reset*" button to restore the factory settings. You can then recalibrate for this magnification.

Once all magnifications have been calibrated, exit this mode by clicking "*Cancel*."


You can retrieve the configured values via "*Management Magnification*":

7.17 Show/Hide Scale Bar

You can find this function under the "Measure" menu. After calibration, the scale can be shown or hidden according to the selected length and unit.


The scale bar will then appear on the display as follows:

7.18 Measure

Under the "*Measure*" menu, there is a wide range of measurement options available. Here is an overview.:

 First, select the image you want to move. The image will now be displayed with a dotted line. Then, mark the cluster of movable points, individual points, or the entire image.

 Select the starting point, click and hold the left mouse button. Move the image to the endpoint, then release the mouse button.

 Select the starting point, click and hold the left mouse button, move the mouse to the endpoint, release the mouse button, and the action is complete (dragging is limited to the horizontal area).

Select the starting point, click and hold the left mouse button, move the mouse to the endpoint, release the mouse button, and the action is complete (dragging is limited to the vertical area).

The first click selects the starting point, the second click defines a straight line, and the third click passes through the third point to create the perpendicular line to the straight line.

Click once to select the starting point, click a second time to define a straight line, and then click a third time to create the perpendicular line to the straight line.

The first click selects the starting point, the second click defines a straight line, and the third click draws both a parallel and a perpendicular line between the two lines.

Click with the mouse on three points at different positions to form an angle.

Select the top-left corner of a rectangle, click and hold the left mouse button, move the mouse to the bottom-right corner of the rectangle, and finally release the mouse button.

Select the center of a circle, click and hold the left mouse button, move the mouse to a point on the arc, and finally release the mouse button.

Click with the mouse on three points at different positions to form a circle.


Select a point with the left mouse button and right-click to finish the operation.

7.19 Measurement settings

Here, you can select settings such as accuracy, font size, and background.

7.20 Export Data


You can export and save the measured data as a CSV file.

	A	B	C	D	E
1	No.	Name	Type	Result	Unit

7.21 Report

You can create a report by selecting "*Export*" in the "*Measure*" menu. The generated DOC file will be exported and can be saved.

The document can then be edited in Word.

S-Viewer Image Analysis Report

Author: _____ Title: _____ Assistant: _____

Sample: _____ Type: _____

Date: 2025-02-10 14:53:58

Measurement:

No.	Name	Type	Result	Unit
1	HL1	HorizontalLine distance	1.000	mm


Comments:

7.22 Annotation

This comment feature allows you to add notes, explanations, or feedback to a specific image. It enables you, as the user, to provide important and additional information.

	Graphics and text can be moved
	Add arrows
	Draw freehand
	Add text
	Add numerical enumeration 1, 2, 3, etc.
	Add X/Y crosshair with full line
	Add X/Y crosshair with a dashed line
	Draw a centerline between two parallel lines
	Add crosshair with subdivisions
	Delete individual measurements, graphics, and text
	Delete all measurements, graphics, and text (with confirmation)
	Cancel measurement

7.23 Save settings

With the "Save Option", you define the parameter information that should be included when saving:

File Name: Prefix + Suffix + File format, e.g., image1.bmp, image2.jpg, etc.

Folder: Select the storage location

Image Type: Choose between Bitmap or JPEG file format


Save with Title: Adds the title from the description into the image

Save with Data/Time: Adds the system date and time into the image


Save with Measure: Writes measurement data into the image

Save with Ruler: Adds the calibration ruler into the image

Description: Includes title and author information in the file metadata


7.24 Language settings

Select a different *language*:

7.25 Arrangement of open windows

Windows in a sequence

Horizontal arrangement of windows

Vertical arrangement of windows

7.26 Additional features

Return to video preview mode

Pause video preview

Capture a single frame and display it in the newly opened window

Interval recording

Video recording

7.27 Interval capture

The 'Interval Capture' dialog box is shown with a blue title bar and a close button. It contains two sections. The top section, titled 'Interval', has three spinners: 'Wait (sec)' set to 0, 'Interval (sec)' set to 1, and 'Images' set to 5. The bottom section has labels for 'Wait (sec)', 'Wait', 'File Name', and 'No.'. At the bottom are three buttons: 'Start', 'Stop', and 'Cancel'.

Wait (sec):

This function will start after the selected number of seconds.


Interval (sec):

Another photo will be taken after the selected number of seconds.

Images:

Number of photos to be taken.

7.28 Video recording

Select path.:

Select the storage location and name for the video files.

Start:

Start video recording.

Pause:

Pause video recording.

Stop:

End video recording; the video will be automatically saved in the selected path.

8 Troubleshooting

8.1 Preview failed, blank screen.

- Check if a light source is present and if the exposure value of the video attribute is normal.
- Check if the camera is connected and the driver is installed in the Device Manager.
- Check the "Input Device" dialog box to ensure that the correct camera is selected.
- Close the software and restart it.
- Close the software, disconnect the camera, reconnect it, and then restart the software.
- Close the software, restart the computer, and then launch the software again.

Please follow the above steps until you can see a normal preview. If the preview still does not appear, please contact your KERN & SOHN service.