

IKA

designed for scientists

The Stirring Bestseller

RW 20 digital and RW 20 digital Package

Now
unbeatable
in Price

%

Reliable. Robust.
Ready for daily lab use.

The RW series is the robust and long-established overhead stirrer line in the market, designed for durability and simplicity with flexible high torque and high-speed operation. This series comes with two gear settings to support high torque or high speed mixing respectively.



WWW.IKA.COM

RW 20 digital

IKA's bestseller overhead stirrer is now available at an even more attractive price, making high-quality stirring technology accessible across industries. The RW 20 digital is a robust and simple overhead stirrer for tough environments and straightforward processes.

YOUR BENEFITS

- › Proven, reliable performance
- › Robust design for demanding environments
- › Easy handling for everyday processes
- › Excellent price-performance ratio



RW 20 digital


Ident. No. 0003593000



Laboratory stirrer designed for **highly viscous applications** and intensive mixing



Extremely powerful laboratory stirrer designed with **high torque**




Digital display for precise monitoring of set and actual speeds

RW 20 digital Package

The RW 20 digital Package is a complete, ready-to-run mixing solution offered at an attractive package price. This high-value, turnkey bundle ensures fast setup, with everything included for immediate use.

Scope of delivery:



RW 20 digital Package

Ident. No. 0009011600

Package
Price
€ 1.188

RW 20 digital

Two speed ranges for universal use from 60 - 2.000 rpm

R 1826

Plate stand

R 1342

Propeller stirrer

R 182

Boss head clamp

Ready to mix? Get your complete RW 20 solution today

Visit www.ika.com or contact your local IKA Sales Team.

IKA-Werke GmbH & Co. KG

Janke & Kunkel-Straße 10, 79219 Staufen, Germany

Phone: +49 7633 831-0, eMail: sales@ika.de



www.ika.com



IKAworlwide // #lookattheblue