



Compatibility Charts

Designed for accuracy,
built for discovery —
empowering life sciences
with precision-driven
equipment

Compatibility Chart -

AHN® myTip® Tecan® Type Tips

Device		Tray Color	Fluent Automation Workstation			Freedom EVO/EVOLyzer		
			LiHa		MCA	LiHa		MCA
			Hanging Tips	SLAS Tips	SLAS Tips	Hanging Tips	SLAS Tips	SLAS Tips
1-304-F50-7	ahn myTip® Robotic Tecan® type 10µL tips Conductive, Flip box	Grey	■	X	X	■	X	X
1-305-F50-7	ahn myTip® Robotic Tecan® type 10µL tips Filter, Conductive, Flip Box		■			■	X	X
1-304-B24-7	ahn myTip® Robotic Tecan® type 10µL tips Conductive, 1-Blister pack		■	X	X	■	X	X
2-304-B24-7	ahn myTip® Robotic Tecan® type 10µL tips Sterile, Conductive, 1-Blister Pack		■	X	X	■	X	X
1-305-B24-7	ahn myTip® Robotic Tecan® type 10µL tips Filter, Conductive, 1-Blister pack		■	X	X	■	X	X
2-305-B24-7	ahn myTip® Robotic Tecan® type 10µL tips Sterile, Filter, Conductive, 1-Blister pack	Green	■	X	X	■	X	X
1-314-F50-7	ahn myTip® Robotic Tecan® type 50µL tips Conductive, Flip box		■	X	X	■	X	X
1-315-F50-7	ahn myTip® Robotic Tecan® type 50µL tips Filter, Conductive, Flip Box		■	X	X	■	X	X
1-314-B24-7	ahn myTip® Robotic Tecan® type 50µL tips Conductive, 1-Blister pack		■	X	X	■	X	X
2-314-B24-7	ahn myTip® Robotic Tecan® type 50µL tips Sterile, Conductive, 1-Blister Pack		■	X	X	■	X	X
1-315-B24-7	ahn myTip® Robotic Tecan® type 50µL tips Filter, Conductive, 1-Blister pack	Blue	■	X	X	■	X	X
2-315-B24-7	ahn myTip® Robotic Tecan® type 50µL tips Sterile, Filter, Conductive, 1-Blister pack		■	X	X	■	X	X
1-329-F50-7	ahn myTip® Robotic Tecan® type 200µL tips Conductive, Flip Box		■	X	X	■	X	X
1-330-F50-7	ahn myTip® Robotic Tecan® type 200µL tips Filter, conductive, Flip Box		■	X	X	■	X	X
1-329-B24-7	ahn myTip® Robotic Tecan® type 200µL tips Conductive, 1-Blister pack		■	X	X	■	X	X
2-329-B24-7	ahn myTip® Robotic Tecan® type 200µL tips Sterile, conductive, 1-Blister Pack	Yellow	■	X	X	■	X	X
1-330-B24-7	ahn myTip® Robotic Tecan® type 200µL tips Filter, conductive, 1-Blister Pack		■	X	X	■	X	X
2-330-B24-7	ahn myTip® Robotic Tecan® type 200µL tips Sterile, Filter, Conductive, 1-Blister Pack		■	X	X	■	X	X
1-404-F50-7	ahn myTip® Robotic Tecan® type 1000µL tips Conductive, Flip box		■	X	X	■	X	X
1-405-F50-7	ahn myTip® Robotic Tecan® type 1000µL tips Filter, Conductive, Flip box		■	X	X	■	X	X
1-404-B24-7	ahn myTip® Robotic Tecan® type 1000µL tips Conductive, 1-Blister pack	Grey	■	X	X	■	X	X
2-404-B24-7	ahn myTip® Robotic Tecan® type 1000µL tips Sterile, conductive, 1-Blister pack		■	X	X	■	X	X
1-405-B24-7	ahn myTip® Robotic Tecan® type 1000µL tips Filter, Conductive, 1-Blister pack		■	X	X	■	X	X
2-405-B24-7	ahn myTip® Robotic Tecan® type 1000µL tips Sterile, Filter, Conductive, 1-Blister pack		■	X	X	■	X	X
1-302-F50-0	ahn myTip® Robotic Tecan® type 10µL tips, clear, Flip Box	Green	■	X	X	■	X	X
1-303-F50-0	ahn myTip® Robotic Tecan® type 10µL tips Filter, clear, Flip Box		■	X	X	■	X	X
1-302-B24-0	ahn myTip® Robotic Tecan® type 10µL tips clear, 1-Blister Pack		■	X	X	■	X	X
2-302-B24-0	ahn myTip® Robotic Tecan® type 10µL tips Sterile, clear, 1-Blister Pack		■	X	X	■	X	X
1-303-B24-0	ahn myTip® Robotic Tecan® type 10µL tips Filter, clear, 1-Blister Pack	Blue	■	X	X	■	X	X
2-303-B24-0	ahn myTip® Robotic Tecan® type 10µL tips Sterile, Filter, clear, 1-Blister Pack		■	X	X	■	X	X
1-312-F50-0	ahn myTip® Robotic Tecan® type 50µL tips clear, Flip Box		■	X	X	■	X	X
1-313-F50-0	ahn myTip® Robotic Tecan® type 50µL tips Filter, clear, Flip Box		■	X	X	■	X	X
1-312-B24-0	ahn myTip® Robotic Tecan® type 50µL tips clear, 1-Blister Pack		■	X	X	■	X	X
2-312-B24-0	ahn myTip® Robotic Tecan® type 50µL tips Sterile, clear, 1-Blister pack	Yellow	■	X	X	■	X	X
1-313-B24-0	ahn myTip® Robotic Tecan® type 50µL tips Filter, clear, 1-Blister Pack		■	X	X	■	X	X
2-313-B24-0	ahn myTip® Robotic Tecan® type 50µL tips Sterile, Filter, clear, 1-Blister pack		■	X	X	■	X	X
1-327-F50-0	ahn myTip® Robotic Tecan® type 200µL tips clear, Flip Box		■	X	X	■	X	X
1-328-F50-0	ahn myTip® Robotic Tecan® type 200µL tips Filter, clear, Flip Box		■	X	X	■	X	X
1-327-B24-0	ahn myTip® Robotic Tecan® type 200µL tips clear, 1-Blister Pack	Grey	■	X	X	■	X	X
2-327-B24-0	ahn myTip® Robotic Tecan® type 200µL tips Sterile, clear, 1-Blister pack		■	X	X	■	X	X
1-328-B24-0	ahn myTip® Robotic Tecan® type 200µL tips Filter, clear, 1-Blister Pack		■	X	X	■	X	X
2-328-B24-0	ahn myTip® Robotic Tecan® type 200µL tips Sterile, Filter, clear, 1-Blister pack		■	X	X	■	X	X
1-402-F50-0	ahn myTip® Robotic Tecan® type 1000µL tips clear, Flip Box	Blue	■	X	X	■	X	X
1-403-F50-0	ahn myTip® Robotic Tecan® type 1000µL tips Filter, clear, Flip Box		■	X	X	■	X	X
1-402-B24-0	ahn myTip® Robotic Tecan® type 1000µL tips clear, 1-Blister Pack		■	X	X	■	X	X
2-402-B24-0	ahn myTip® Robotic Tecan® type 1000µL tips Sterile, clear, 1-Blister pack		■	X	X	■	X	X
1-403-B24-0	ahn myTip® Robotic Tecan® type 1000µL tips Filter, clear, 1-Blister pack		■	X	X	■	X	X
2-403-B24-0	ahn myTip® Robotic Tecan® type 1000µL tips Sterile, Filter, clear, 1-Blister pack	Yellow	■	X	X	■	X	X
	pure nested conductive tips 10ul		X	■	X	X	■	X
	pure nested conductive tips 10ul with filter		X	■	X	X	■	X
	sterile nested conductive tips 10ul with filter		X	■	X	X	■	X

■ compatible X not compatible

Updated as on Oct 2025. For latest compatibility chart scan here

The compatibility chart above is for your reference. We are not responsible for technical modifications or changes made to the mentioned PCR devices.



Compatibility Chart -

AHN[®] myTip[®] Tecan[®] Type Tips

LiHa		MCA	Device
Hanging Tips	SLAS Tips	SLAS	
■	x	x	Abbot Diagnostics M2000
■	-	-	Adaltis NexGEN Series
■	■	■	Aurora Biomed VERSA
■	x	x	AusDiagnostics Ultrplex Alliance system
■	x	x	AXA Diagnostics Everex - SkyLAB 752 -Italy
■	x	x	Corbett CAS-1200
■	x	x	Cybio Analytik Jena CyBio FELIX
■	x	x	DaAN- China Gene Automated Nucleic Acid Extraction system
■	x	x	Lonza- Swiss Pyrotec PRO Automated Robotic system
■	x	x	Hologic - USA Panther Fusion
■	x	x	Mobidiag Hologic Group - Finland Amplidiag Easy system
■	■	x	Perkin Elmer Evolutiontm P3
■	■	x	Janus
■	x	x	Qiagen QIAgility
■	-	x	Grifols Procleix Panther system
-	■	-	Yourgene Health Yourgene® Sp150
■	-	-	ThermoFiser Scientific Thermo Fisher Scientific Amplitude Solution
■	■	-	Sansura Biotech Sansura Biotech Natch S
■	-	-	Mobidiag Mobidiag Amplidiag® Easy
-	■	-	Miltenyi Biotec Miltenyi MultiMACS™ X
■	-	-	NGNY Devices NGNY AQUA7000
■	-	-	TelesisBio BioXp 3200/3250/9600

■ compatible x not compatible - Under Validation

Updated as on Oct 2025. For latest compatibility chart scan here

The compatibility chart above is for your reference. We are not responsible for technical modifications or changes made to the mentioned PCR devices.



Distributed by:



ahn[®]

AHN Biotechnologie GmbH
Uthleber Weg 14
99734 Nordhausen Germany
P: +49(0)3631/65242-0
F: +49(0)3631/65242-90
E: contact@ahn-bio.com

www.ahn-bio.com

