



Compatibility Charts

Designed for accuracy,
built for discovery —
empowering life sciences
with precision-driven
equipment

Compatibility Chart - Tips/Filter Tips

■ compatible x not compatible

Updated as on Oct 2025. For latest compatibility chart scan here

The compatibility chart above is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes. Specifications like inaccuracy and imprecision values have not been tested.

[illegible]

Compatibility Chart - Tips/Filter Tips

■ compatible x not compatible

Updated as on Oct 2025. For latest compatibility chart scan here

The compatibility chart above is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes. Specifications like inaccuracy and imprecision values have not been tested.

[illegible]

Compatibility Chart - Tips/Filter Tips

Compatibility Chart - Tips/Filter Tips						
Sartorius	HTL (Coming)	Hamilton	Nichiryo	Socorex	IKA	Primo
	Labmate Pro (some as Labnet BioFette)/ Pipette DV Series manual single and multi-channel pipettes	Softgrip manual single channel and multi-channel pipette	Nichipor Premium / EX II / EX I Multi manual single and multi-channel pipette	Socorex Acura manual single and multi channel pipette	IKA PETTE Vario / Multi / Fix manual single and channel pipette	Primo Pipettes single channel and multi channel pipette
	(Proline/Milne/Tacta) manual single channel (Milnes3/Teille/acta) Manual Multi channel (electronic single channel (eIine/Picus)) electronic multi channel (Picus)	manual multi	Manual Single channel	multi channel	single channel	single channel
Cat. No.	Description	AHN myTip® Pipette Tips				
1-002-81-0 AHN myTip® 10 µL clear						
1-002-C5-0 AHN myTip® 10 µL clear						
1-003-C5-0 AHN myTip® 10 µL clear, sterile						
1-062-81-0 AHN myTip® 10/20 µL XL Clear						
1-062-C5-0 AHN myTip® 10/20 µL XL Clear						
1-063-C5-0 AHN myTip® 10/20 µL XL, clear, sterile						
1-121-86-0 AHN myTip® 200 µL clear						
1-121-C5-0 AHN myTip® 200 µL clear						
1-122-C5-0 AHN myTip® 200 µL clear, sterile						
1-121-86-1 AHN myTip® 200 µL yellow						
1-121-CS-1 AHN myTip® 200 µL yellow						
1-122-C5-1 AHN myTip® 200 µL yellow, sterile						
1-130-70-0 AHN myTip® 200 µL Gelloader, clear						
1-130-C4-0 AHN myTip® 200 µL Gelloader, clear						
1-131-C4-0 AHN myTip® 200 µL Gelloader, clear, sterile						
1-127-80-0 AHN myTip® 300 µL clear						
1-127-CS-0 AHN myTip® 300 µL clear						
1-128-CS-0 AHN myTip® 300 µL clear, sterile						
1-209-80-2 AHN myTip® 1000 µL Short, clear Blue						
1-202-80-0 AHN myTip® 1000 µL clear						
1-202-C4-0 AHN myTip® 1000 µL clear						
1-203-C4-0 AHN myTip® 1000 µL clear, sterile						
1-202-80-2 AHN myTip® 1000 µL blue						
1-202-C4-2 AHN myTip® 1000 µL blue						
1-203-C4-2 AHN myTip® 1000 µL blue, sterile						
1-204-80-0 AHN myTip® 1250 µL clear	X					
1-204-C4-0 AHN myTip® 1250 µL clear	X					
1-205-C4-0 AHN myTip® 1250 µL clear, sterile	X					
AHN myTip® Mat Macro Pipette Tips						
1-218-61-0 AHN myTip® Mat 5 mL clear			X		X	X
1-222-61-0 AHN myTip® Mat 5 mL clear			X		X	X
1-210-84-0 AHN myTip® Mat 5 mL clear			X		X	X
1-212-84-0 AHN myTip® Mat 5.5 mL clear			X		X	X
1-216-84-0 AHN myTip® Mat 5.5 mL clear, sterile			X		X	X
1-213-50-0 AHN myTip® Mat 10 mL clear				X	X	X
1-220-50-0 AHN myTip® Mat 10 mL clear			X		X	X
1-224-50-0 AHN myTip® Mat 10 mL clear, sterile			X		X	X

■ compatible x not compatible

Updated as on Oct 2025. For latest compatibility chart scan [here](#)



The compatibility chart above is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes. Specifications like inaccuracy and imprecision values have not been tested.

Compatibility Chart - Tips/Filter Tips

Compatibility Chart - Tips/Filter Tips											
Description	Cat. No.	Sartorius				HTL (Coming) Labmate Pro (same as Labnet BioPipette) / Pipette DV Series	Hamilton Softgrip manual single and multi channel pipette	Nichiyo Nichiyo Premium / EX II / EXII Multi manual single and multi channel pipette	Socorex Socorex Acura manual single and multi channel pipette	IKA IKA PETTE Vario / Multi / Fix manual single and multi channel pipette	Primo Primo Pipettes single channel and multi channel pipette
		Proline Plus / Picus / Biohit / MLine / Tacta / e300 manual single channel and electronic single and multi channel pipettes									
		Manual Multi (Mline3/Proline/T acta)	electronic single channel (eline/Picus)	electronic multi channel (Picus)							
AHN myTip® Filter Tips											
AHN myTip® FT 10 µL, clear	2-002-81-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 10 µL, clear, sterile	2-003-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 10/20 µL, clear	2-062-81-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 10/20 µL, clear, sterile	2-063-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 20 µL, clear	2-117-80-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 20 µL, clear, sterile	2-118-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 100 µL, clear	2-119-80-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 100 µL, clear, sterile	2-120-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 200 µL, clear	2-127-80-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 200 µL, clear, sterile	2-128-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 300 µL, clear	2-134-80-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 300 µL, clear, sterile	2-135-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 1000 µL, clear	2-202-80-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 1000 µL, clear, sterile	2-203-C4-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 1250 µL, clear	2-204-80-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® FT 1250 µL, clear, sterile	2-205-C4-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT Low Retention Pipette Tips											
AHN myTip® LT 10 µL, clear	4-002-81-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT 10 µL, clear, sterile	4-003-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT 10/20 µL, clear	4-062-81-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT 10/20 µL, clear, sterile	4-063-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT 200 µL, clear	4-121-86-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT 200 µL, clear, sterile	4-122-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT 300 µL, clear	4-127-80-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT 300 µL, clear, sterile	4-128-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT 1000 µL, clear	4-202-80-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT 1000 µL, clear, sterile	4-203-C4-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT 1250 µL, clear	4-204-80-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LT 1250 µL, clear, sterile	4-205-C4-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LFT Low Retention Filter Tips											
AHN myTip® LFT 10 µL, clear, sterile	5-003-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LFT 20 µL, clear, sterile	5-063-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LFT 20 µL, clear, sterile	5-118-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LFT 100 µL, clear, sterile	5-120-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LFT 200 µL, clear, sterile	5-128-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LFT 300 µL, clear, sterile	5-135-C5-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LFT 1000 µL, clear, sterile	5-203-C4-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel
AHN myTip® LFT 1250 µL, clear, sterile	5-205-C4-0					0.5-10 µL	0.2-2 µL	1-10 µL	5-50 µL	30-300 µL	multi channel

compatible not compatible

Updated as on Oct 2025. For latest compatibility chart scan here



The compatibility chart above is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes. Specifications like inaccuracy and imprecision values have not been tested.

Distributed by:



ahn[®]

AHN Biotechnologie GmbH
Uthleber Weg 14
99734 Nordhausen Germany
P: +49(0)3631/65242-0
F: +49(0)3631/65242-90
E: contact@ahn-bio.com

www.ahn-bio.com

