

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR® Culture Tubes With Screw Thread
Country of Origin	Manufactured in India
Date of Issue (yyyy-mm-dd)	2024-07-23

Catalogue Number	Thread	Height [mm]	Diameter [mm]
212-0849	13-415	100	13
212-0850	15-415	100	16
212-0851	15-415	125	16
212-0852	15-415	150	16
212-0853	18-415	125	20
212-0854	18-415	150	20

Screw caps for culture tubes.

Catalogue Number	Description	Thread	Colour
212-0855	Screw cap with bonded rubber liner	13-415	Black
212-0856	Screw cap with seal liner	15-415	Black
212-0857	Screw cap with bonded rubber liner	18-415	Black

Quality System Compliance

Products are manufactured under the **ISO 9001:2015 & ISO 18001:2007** standard. Products are inspected and controlled through whole production processing in accordance with current applicable product specifications and QC SOP.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

The product is in conformance with the requirements of the following standards and regulations

Technical standards for products: NA

Calibration procedure: NA

Material	Culture Tube: Borosilicate Glass Screw Cap: Phenolic, Natural Rubber.
Sterilization	Not applicable.
ATP Assay	Not Applicable.
DNase & RNase Free	Not Applicable.
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable.
Non-Pyrogenic Statement	Not Applicable.
Latex Statement	Latex is used as raw material or in manufacture process (which is not applicable for culture tubes).
BPA Statement	Not Applicable.
DEHP Statement	Not Applicable.
RoHS	Not Applicable.
REACH Statement	Product meet the criteria for authorization under the European Union's REACH regulation (EC) No. 1907/2006. To the best of our knowledge product does not contain SVHC's above the regulatory threshold of 0.1%.
Storage Conditions	Store at room temperature in normal warehouse conditions.
Shelf Life	No predefined shelf life.

Disclaimer: VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. The purpose of the product is for use in laboratory only.
This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoQ using new layout	07-2024