

Certificate of Conformity

VWR Catalogue Number	See List
Description	VWR®, Centrifuge Tubes, PTFE
Country of Origin	Manufactured in United Kingdom
Date of Issue (yyyy-mm-dd)	2025-11-14

Round Bottom

Catalogue Number	Capacity [mL]	Diamention Ø×L [mm]	Closure Type
525-0178	8	12×100	Without cap
525-0179	13	16×100	Without cap
525-0180	17	18×100	Without cap
525-0181	50	25×145	Without cap
525-0182	85	40×100	Without cap

Quality System Compliance

This certificate confirms that the product has been sourced from a supplier or manufacturer who has declared their Quality Management System (QMS) complies with **ISO 9001:2015**.

QC Testing

Representative products samples are collected and inspected in accordance with current applicable QC requirements.

Product Specifications

Material	Virgin Polytetrafluoroethylene (PTFE). Raw Material complies with Directive FDA 21 CFR 177.1550 and complies with USP Class VI.
Temperature Resistance	–200 °C to +280 °C.
Sterilization	Non-Sterile.
ATP Assay	Not Applicable

DNase & RNase Free	Not Applicable
BSE/TSE	No use of any raw material produced from or substances derived from animal origin. The manufacturing process of the product does not use any ingredient of animal origin and no material derived from or exposed to animals affected by or under quarantine for transmitting Animal Bovine Spongiform Encephalopathy/Spongiform Encephalopathy (BSE)/(TSE).
Cytotoxicity	Not Applicable
Non-Pyrogenic Statement	Not Applicable
Latex Statement	This product is not made from natural rubber latex.
BPA Statement	Bisphenol A absent Bisphenol not used in the manufacture of the raw material and are not expected to be present.
DEHP Statement	The substance Phthalates are not used or intentionally incorporated in the manufacture of the raw materials or final product.
RoHS	No substances (Lead, Cadmium, Mercury, Hexavalent Chromium (Cr6+), Poly Brominated Biphenyls (PBB), Poly Brominated Diphenyl ethers (PBDE), Benzyl Butyl Phthalate (BBP), Dibutyl Phthalate (DBP), Diisobutyl Phthalate (DIBP)) are used in manufacturing the raw materials and final product. No routinely analyse is performed.
REACH Statement	No substances as listed on >Substance of Very High Concern< candidate list are used in manufacturing of the raw materials used to produce this product. No routinely analysis is performed.
Storage Conditions	Ambient Temperature.
Shelf Life	Indefinite when stored under ideal conditions.

Disclaimer: VWR's Quality Management System (QMS) is ISO 9001:2015 certified as a wholesale distributor of finished products, ensuring effective oversight of finished product handling, traceability, and distribution processes. VWR is not a manufacture and does not engage in manufacturing, co-manufacturing or processing activities such as design, packaging, labeling, testing or quality release. Consequently, the supplier or manufacturer's QMS certification may differ from VWR's.

VWR states that this declaration will not discharge the user from their obligation to ensure the product is suitable for the intended use. This product is intended for research use only. It is not for use in diagnostic procedures, therapeutic applications, or any human or clinical use.

This document has been produced electronically and is valid without a signature.

Version History

Sl.no	Version number	Description	Month & Year
1	1.0	Created General CoC using new layout	2025-11