ahn myVial®

Diagnostic Vial



Fully Automated Manufacturing

Manufactured from 100K class cleanroom Constructed from USP Class VI US FDA Approved PETG material

Exceptional Clarity

Ultra-clear PETG ensures effortless visual inspection of contents, making it ideal for diagnostic and clinical applications.

Diverse Volume Ranges

Available in three volumes 5ml, 10ml, 20ml



Leakproof Security

Each vial is fitted with a precision-engineered HDPE screw cap and an integrated polyethylene cap liner that ensures a secure seal, minimizing leakage during transport.

Non-Cytotoxic & Non-Pyrogenic

Tested and certified to be safe for biological samples, supporting sensitive diagnostic and pharmaceutical workflows.

Sterility Assurance

Gamma sterilized (SAL 10^{-6}) in ISO Class 7 cleanroom conditions for assured contamination-free performance.

Certified

RNase, DNase, Endotoxin, PCR Inhibitor-free

Packaging

Available in a variety of standard sizes with tailored packaging options to match your workflow needs

AHN Biotechnologie GmbH
Uthleber Weg 14
99734 Nordhausen
Germany
www.ahn-bio.com

ahn

P: +49(0)3631/65242-0
F: +49(0)3631/65242-90
E: contact@ahn-bio.com

ahn myVial®

Diagnostic Vial

Ordering information

Cat Code	Decription	Quantity/Case
AHN-PETG-10-S	ahn® myVial Diagnostic Vial PETG 10 ml, sterile, racked	5 x 20 pcs
AHN-PETG-10-S-B	ahn® myVial Diagnostic Vial PETG 10 ml, sterile, racked, bulk	5 x 100 pcs
AHN-PETG-20-S	ahn® myVial Diagnostic Vial PETG 20 ml, sterile, racked	5 x 20 pcs
AHN-PETG-20-S-B	ahn® myVial Diagnostic Vial PETG 20 ml, sterile, racked, bulk	5 x 100 pcs
AHN-PETG-5-S	ahn® myVial Diagnostic Vial PETG 5 ml, sterile, racked	5 x 20 pcs
AHN-PETG-5-S-B	ahn® myVial Diagnostic Vial PETG 5 ml, sterile, racked, bulk	5 x 100 pcs

AHN Biotechnologie GmbH
Uthleber Weg 14
99734 Nordhausen
Germany
www.ahn-bio.com

ahn

P: +49(0)3631/65242-0
F: +49(0)3631/65242-90
E: contact@ahn-bio.com