ahn *reagent* Reservoirs

Construction - formulated from of a unique blend of polystyrene

Extra large, sturdy base with broad edges - enables it to stand fir mLy on your work bench

Graduation marks - helps identify any remaining liquid and facilitates easy, accurate dispensing

Angled interior - helps the collection of the liquid and decreases reagent waste

Spillproof - no need to worry about sample spillage while you are processing liquids

Easy pouring - to store any remaining liquid you can empty the reservoir using the corner spouts

Compact storage - easy to stack and the lip at the bottom of the reservoir allows them to be easily separated from each other



Ideal for using with all AHN Single and Multi-channel Pipettes

- our reservoirs come in three different sizes, to suit various pipettes:

1. 10 mL, V-shaped bottom (ideal with 8- and 16-channel pipettes)
2. 25 mL, flat bottom (ideal with multi-channel pipettes)
3. 50 mL, V-shaped bottom (ideal with 12-channel pipettes)



Ordering information

Description	Colour	Pack Type	Sales Unit	Cat. No.
Reagent Reservoirs **10 mL** 	white	bag	case / 5x100	3-701-11-9
Reagent Reservoirs **25 mL** 	white	bag	case / 5x100	3-702-11-9
Reagent Reservoirs **50 mL** 	white	bag	case / 5x100	3-703-11-9

+49(0)3631/65242-0 contact@ahn-bio.com www.ahn-bio.com