ahn *my*Tube®


Spin Column System
SC

Spin column system - consists of a 0.8 mL filter tube (round bottom design) equipped with corresponding filter (see filter options below) and a 2.0 mL receiver tube

Elimination of risk of contamination or loss of sample - high-quality raw material meets purity standards, highly resistant to extreme mechanical strain and has excellent chemical resistance

Centrifuge specifications: spins up to 25000 x g with a fixed angle rotor, or up to 30000 x g with a swing out rotor





ahn *my*Plate®

Microfilterplate System
FP


Microfilterplate - 8x12 matrix, 96 individual drip directors, 1.0 mL capacity, in accordance with SBS-Standard. Equipped with the corresponding filter (see filter options below)

Elimination of any risk of sample contamination - manufactured from absolutely clear high-quality virgin polystyrene

Filter Options	SC	FP
Glass fibre (GF) filter characteristics - recommended for the purification of DNA and RNA / prefiltering of dirty solutions / rapid flow rate and high particle loading capacity	■	■
Polyethylene (PE) filter characteristics - suitable for particle-removing coarse filtration / can be used both as a prefilter or as support filter in combination with other filters / hydrophobic / very good chemical resistance	■	■

+49(0)3631/65242-0 contact@ahn-bio.com www.ahn-bio.com

Ordering information



Description - spin column system		Colour	Pack Type	Sales Unit	Cat. No.
myTube® SC **0.8 mL** filter tube, **GF-F2 filter, 2.0 mL** receiver tube	1	clear	bag	2x100	3-320-20-0
myTube® SC **0.8 mL** filter tube, **GF-N2 filter, 2.0 mL** receiver tube	1	clear	bag	2x100	3-321-20-0
myTube® SC **0.8 mL** filter tube, **PE 20-60 µm filter, 2.0 mL** receiver tube	1	clear	bag	2x100	3-340-10-0
Description - individual components spin column system		**Colour**	**Pack Type**	**Sales Unit**	**Cat. No.**
myTube® SC **0.8 mL** filter tube w/o filter, w/filter ring	2	clear	bag	2x1000	3-331-50-0
myTube® RT **2.0 mL** receiver tube	3	clear	bag	10x1000	3-202-80-0
Description - microplates		**Colour**	**Pack Type**	**Sales Unit**	**Cat. No.**
myPlate® FP 96-wells, **1.0 mL, GF-FF** filter	4	clear	bag	5x5	3-210-25-0
myPlate® FP 96-wells, **1.0 mL, GF/N 1.6 µm** filter	4	clear	bag	5x5	3-210-25-2
myPlate® FP 96-wells, **1.0 mL, PE 20-60 µm** filter	4	clear	bag	5x5	3-225-25-0

Further filter options are available on request:

- Nylon - recommended for filtering protein-free culture media
- Cellulose Acetate - recommended for applications requiring low protein binding
- Regenerated Cellulose - suitable for particle removal from solvents
- Polyvinylidene Fluoride - used in general filtration and sample preparation applications

Please contact contact@ahn-bio.com